

MID3: Ligation & Transformation

9/24/13

1. Pre-lab discussion
2. Ligation -- scale as needed
50 ng bKb
3. Get rid of the salt!
4. Leslie here to talk about Abstracts
5. Transform e.coli and plate on LB/Agar + Ab

Review MID2 FNT: Methods section

going forward. - PCR (use your primers!)

- Sub-cloning in pCX-WNX \Rightarrow pCX-4326
FO

① What does the reader need to know?

- primer sequences - text cues

- PCR cycling conditions

- [final] * No μ h



- RE: U(nits) of activity - HF EOR1
- XbaI

② Technical language: - flap/landing
- "clean" \rightarrow purification
- "run on a gel" \rightarrow evaluated.
- separated

③-4 Sub-section
(more specific)

* Introductory sentence

Step I: Build the system!



$$6 \text{ kb} \rightarrow \frac{84 \text{ ng}}{20 \mu\text{L}} = \frac{4.2 \text{ ng}}{\mu\text{L}} \times 4 \approx 16 \frac{\text{ng}}{\mu\text{L}} \approx 3 \mu\text{L } 6 \text{ kb}$$

How to estimate the correct volumes for the ligation:


50 ng 6 kb

$$500 \frac{\text{Da}}{\text{bp}} \times 2 = \left(1000 \frac{\text{Da}}{\text{bp}} \right) (4200 \text{ bp}) \approx 4.2 \times 10^6 \text{ g/mol}$$

50 ng	nmol	$\approx 12 \text{ fmol} \times 4 \approx 4.8 \text{ fmol insert}$ $\Downarrow \text{ ng} \Rightarrow \mu\text{L}$
	$4.2 \times 10^6 \text{ ng}$	

Overview: Ligation


What effects the efficiency of ligation?



• [salt] — 50mM Tris
Mg²⁺, ATP, DTT
↑
availability of ATP

• temperature
16°C O/N
25°C (RT) 10 min

Your Ligations

	bkb + insert, no ligase	[bkb only, + ligase]	bkb + insert, + ligase
What does this control for?	bkb uncut	Single cut bkb	Expt 
pCX-NNX bkb	? uL	? uL	? uL ^{50ng}
Δ5 PCR product	? uL	xxx	? uL
10x buffer	1.5 uL	1.5 uL	1.5 uL
T4 DNA Liagase	xxx	0.5 uL	0.5 uL
Water*	to 15 uL	to 15 uL	to 15 uL

*not including enzyme volume

XL-1 Blue

Overview: Transformation



chemically competent


CaCl₂




ice

heat shock @ 42°C - 90sec

LB/Agar + Amp




Left Pit



Right Pit


division

clonal growth



37°C
LB

(No Amp)



harvest DNA
MDS

<http://www.yourwildlife.org/2012/08/not-all-pits-are-equal/>

Your Transformations

Tube	Transformation	Expectation:	What if?
Not doing	nothing (just the plate)	∅	lots? - contamination
1	pCX-EGFP ^(5ng)	A lot 100's	∅? - overheated - weren't competent - too Ab
2	bkb + insert; no ligase	∅ - 1	A lot? - contaminating uncut plasmid
3	bkb; + ligase	10's - 100	lots? 11 ^ single cut
4	bkb + insert; + ligase	50 - 100 > 200 - 400	Nothing? - transformation - <u>ligation</u>

transformation efficiency:

$\frac{\# \text{ of colonies}}{\text{amt of DNA}}$

(Not plating everything)

Today in the lab:

- Set up ligations using your calculations from the FNT -- remember that total volume of bkb + insert cannot be greater than 13.5 uL
- Clean up ligation -- talk about Abstracts
- Transform into e.coli and then plate (with fire)

Next time in the lab:

- Minipreps
- Diagnostic digests
- Intro to Tissue Culture!

