

Now the fun begins!

MIDI: Microbial DNA Extraction

2/6/13

Announcements

- A few notes about the lab practical
- Please sign up for EN and share notebook by end of today

Lab Notebooks

Module Overview

Lab Overview

FNT

Lab Notebooks: Evernote (EN)

• Known issues
- browser version ≠
stand alone version
(Share/save)

The screenshot displays the Evernote web interface. At the top, the Evernote logo and a search bar are visible. The user's name, 'shannon_hughes_20109', is shown in the top right corner. The main content area is divided into three sections: a left sidebar with 'Shortcuts', 'Notebooks', and 'Tags'; a central 'All Notes' list; and a right-hand editing area for the selected note. The note, dated 'Wednesday, 9/4/13', contains a bulleted list of topics: 'Orientation lecture -- last four slides', 'Pre-lab lecture -- add EN and 'new stuff' slide', 'Demo of Evernote', and 'Better explanation of new safety workflow'. The interface includes a rich text editor with various formatting options like font size, bold, italic, and underline. A 'Download Now' button for the Evernote Web Clipper is located in the bottom left corner of the interface.

From protocol to lab notebook

1. Begin by adding the correct amount of water to a 200 ul PCR tube. Add that amount +1 ul to a second PCR tube.
2. Next add the primers to each reaction. Be sure to change tips between additions.
3. Next add template to the first reaction tube.
4. Finally add PCR Master Mix to each tube, pipetting up and down to mix. Leave your tubes on ice until the entire class

Statement of purpose: Today we will design primers to delete 32 bp from the 5' end of GFP and flank the sequence with new restriction sites. Then we will prepare truncated GFP by PCR as an insert for later cloning.

Design primers for GFP insert (M1D1 Part 1)

See attached Word document.

PCR to make GFP insert (M1D1 Part 2)

Added 27 uL H₂O to expt'l, 28 uL H₂O to control sample.


Added [1 uL] primer and [20 uL] Master Mix (last) to both samples, and 1 uL template to expt'l only!

Rxn ready at 3 pm → on ice → thermal cycler started at 4 pm.

Thanks to Agi Stachowiak for this slide!

Module I Conceptual Overview

Mod 1 \Rightarrow 2 independent experiments



Real World Context -- Bird Microbial Communities

What is our primary research question?

How does the microbiome (bacteria population) differ/vary w/ environment?

- bird type
- male/female
- location

What are the broader impacts of our research?

What accounts for
flu susceptibility?

You will amplify the 16S rRNA gene to profile the microbiome of New England gulls

two types of gulls → ~~Ring-billed~~ Ring-billed Herring gull

Sources:



South Bay Center parking lot - South Boston
Jan 4, 2014 (Photo by Darren McColester/Getty Images)



Carson Beach, South Boston (1940): John Sanroma -
Boston Public Library Flickr Stream

Samples:

cloacal samples




^ cloacal
KISS

Male/
Female



Undisclosed to preserve bird privacy.

Bird Microbial Communities -- Experimental Overview



Today: Purify DNA from coacal samples

Both → Red, yellow, green, pink

Different → orange, blue


★ cross contamination (minimize) ★
- filter tips
- tube changes

Inhibitors of Polymerase Chain Reaction (PCR) in coacal samples:

Inhibit enzyme activity {

- heme - bilirubin
- bile salts
- Humic cpds → degradation products


DNA Extraction from Stool: Part I



DNA Extraction from Stool: Part III



Qiaprep column:
Silica Resin



[Promega.com]

- 1) chaotropic salt [high]
(broken H⁺ bonds)
- 2) bind silica - ethanol
- 3) elute - "high" pH ~ 8.0
- low salt


Stool Lysate




Bind




Wash



Elute



Purified DNA



Today in Lab:

Waste disposal: save all tubes, rinse 2-3x with water wash bottle over marked waste stream in fume hood (safety glasses!).

Step 1: Stool lysis through adding enzymes

- Keep track of tube changes ~45 min
- Use filtered pipette tips — share.

Step 2: 1.5 hr incubation

- *Lab practical*
- Prepare tubes for later steps

Step 3: DNA purification using silica resin column

FNT: wiki page, MID3 paper