


MID3: PCR and Paper Discussion


2/15/13

Announcements

- First lab treat:



Another source of DNA
for microbiome studies



- FNT -- long due to Holiday
- PCR Review and Gel Electrophoresis
- Set up PCR
- Atissa will be here, then journal club!

FNT Assignment

1. Larger assignment than usual. Changed a bit since Wednesday.
 - Draft methods section on MID2-MID3
 - Partial introduction for microbiome work
 - Primer design table + thoughts
 - Note: Methods & Intro due on Stellar
2. I will be out of town this weekend -- contact me through gmail account
3. Office hours next week: Monday, 8pm & Tuesday, 11am

HOMEWORK →

General

[edit topic](#) - [delete topic](#) - [add a](#)

[Due MID4, T/R](#) [edit](#) - [delete](#)
Due 21 February 2013 1:05 p.m.

[Due MID4, W/F](#) [edit](#) - [delete](#)
Due 22 February 2013 1:05 p.m.

[Due MID5, T/R](#) [edit](#) - [delete](#)
Due 26 February 2013 1:05 p.m.

[Due MID5, W/F](#) [edit](#) - [delete](#)
Due 27 February 2013 1:05 p.m.

[Due MID7, T/R](#) [edit](#) - [delete](#)
Due 05 March 2013 1:05 p.m. Pc

[MID7, W/F](#) [edit](#) - [delete](#)
Due 06 March 2013 1:05 p.m. Pc

Methods Section Tips


- Divide into sub-sections!
 - Put in a logical order -- list primer THEN list PCR reagents, etc
 - Start with a overview sentence “16S rRNA genes were amplified from each sample...”
- Methods are clear and concise explanations
 - The methods section is not a benchtop protocol**
 - Space-wise, avoid tables/lists when a sentence will do
 - Sentence-wise, avoid extra words
 - Content-wise, cover what's needed and only that needed to understand and replicate your work
- Think about the most flexible units -- concentration vs volume. (For replication!)

Methods section exercise

- Consider the following passage: “Template DNA (5 ng) and primers were mixed with 20 uL of 2.5X Master Mix in a PCR tube. Water was added to 50 uL. A tube without template was prepared and labeled control.”
 - What information is missing?
 - What information can be cut?

Thanks to Agi for this slide!

Bird Microbial Communities -- Experimental Overview



PCR


Component	Function
template DNA	Original Copy
polymerase Tag ← Pfu	Catalyzes DNA addition
dNTPs	Building Blocks
primers	★ Select and initiate new sequence
Buffer Mg ²⁺	Optimal chemical environment

BSA


★ Hot Start

DNA Electrophoresis (EP): Principle

Agarose gel



DNA



Agarose and DNA are both *polymers*

Driving force for separation:


DNA moves $-$ to $+$ because of *phosphate groups*

Separation is according to: *Size*

Smaller

DNA moves faster because


porosity



$+$

Thanks to Agi for this slide!


How do we visualize the DNA?



Loading dye

- glycerol
- C
- xylene cyanol
- RNase

← Sybr green



sr.wikipedia.org




awkward photo with your advisor

DNA EP: Analysis

MW


IBS

N.T.C



DNA ladder:

Relationship:



From MIT OCW

About *V. corneae*... and your primers

- Gave you poor reference sequence
 - Broad = shotgun sequence, avoided it for simplicity
 - DaSilva complete: 3 mismatch, 1 gap
 - Baker complete: 2 mismatch, 25(!) gaps
- Gave you wrong direction
 - Target T_a , not T_m , of 58 °C is best (matches V1/PMP2)
- What we will do
 - revise your primers to Broad seq., checked new T_m & G/C
- Your options now
 - accept our revision (we will run a second PCR at 53 °C)
 - revise your primers to a T_m of 63 °C

Thanks to Agi for this slide!

Today in Lab (M1D3)

- Set up PCR rxns
 - Change pipet tips between samples, primers, etc.
 - Keep PCR tubes cold!
 - Write small *directly* on the PCR tubes – do not put sticky labels in the PCR machine.
- Discuss paper from writing POV ~2 pm
- Presentation on giving talks from Atissa ~2:20
- Polish your slide ~ 3:15
- Discuss paper from technical POV *and* get feedback about your slide ~ 3:30-5 pm

Thanks to Agi for this slide!