


# M1D3: Measure DNA damage with CometChip

09/20/16

1. Quiz
2. Pre-lab Discussion
3. CometChip loading discussion
4. Load CometChip, Induce DNA damage
5. CometChip electrophoresis

# Overview of


## “M1: Measuring Genomic Instability”



### 1. Optimize comet chip assay

- Test loading variables


discuss results today



### 2. Use comet chip assay to measure DNA damage / repair

- Measure effects of MMS and  $H_2O_2$  on BER
- Assess repair variability in healthy individuals →


Next week



### 3. Use immuno-fluorescence assay to visualize DNA repair

- Examine effect of  $H_2O_2$  on DSB abundance


# Start with CometChip from M1D1...




# Treat cells with 6 different doses of H<sub>2</sub>O<sub>2</sub>



Concentrated H<sub>2</sub>O<sub>2</sub> is a mutagen!  
Wear gloves and a lab coat at all times.



← H<sub>2</sub>O<sub>2</sub> treatment  
(in macrowells of  
96-well plate)



} 3 replicates  
← 6 doses

$$C_1 V_1 = C_2 V_2$$


## Treat cells with 6 different doses of H<sub>2</sub>O<sub>2</sub>

- Stock H<sub>2</sub>O<sub>2</sub> is 10M
- Add 100μl of drug dose to each macrowell
- Add drug to the “live” chip and “dead” chip in parallel
- Prepare dilutions in cold PBS
- Treat with : 0, 10, 20, 40, 60, 80 μM
- *Dilutions should minimize waste!*

1) 10mL PBS + 10μl Stock H<sub>2</sub>O<sub>2</sub> = 10 mM

2) 10mL PBS + 100μl ① = 100 μM

make 1ml  
(1000μl) of  
each dose





# Why analyze dead cells?

enzyme dead

live

dead

DNA damage




no tail=  
no damage

tail=  
DNA damage

Note: W/F students will treat their TK6 cells with another drug, MMS.  
➤ Will the two chemicals have the same DNA damaging effect?


H2O2 vs MMS



# How do $\text{H}_2\text{O}_2$ and MMS damage DNA?

oxidative damage


$\text{H}_2\text{O}_2$



Repaired  
by  
BER

alkylation damage

MMS



# Lysis & electrophoresis of CometChips

- Alkaline lysis solution (> 1 h at 4°C)
  - 2.5 M NaCl, 100 mM Na<sub>2</sub>EDTA, 10 mM Tris, (water)
  - pH 10
  - Triton X-100 *detergent*

} breaks open cells  
denatures proteins

- Electrophoresis buffer (40 min incubation + 30 min run, at 4°C)
  - 0.3 M NaOH, 1 mM Na<sub>2</sub>EDTA
  - pH 13.5 *very high pH*
  - gellbond adhered with double-sided tape


} unwinds  
DNA

- Electrophoresis
  - 30 min
  - 1 V/cm (or 16 V for us)
  - 300 A

} very slow

- Neutralize & dye
  - Tris Base solution ~ pH 7.4 (biological)
  - SYBR Gold

→ DNA stain





# Tips to write Methods (due M1D4)

- Include enough information to replicate the experiment
  - list manufacturers information
  - Be **concise and clear** in your description
- Use subsections with descriptive titles
  - Put in logical order
  - Begin with topic sentence to introduce purpose
- Use clear and concise full sentences
  - Eliminate tables and lists
  - Passive voice expected
- Use the most flexible units
  - Write concentrations (when known) rather than volumes
- Eliminate '109 specific details
  - Do not include details about tubes and water!
  - Assume reader has some biology experience

# Methods Practice

No 1 correct way to write a methods section!!

Maintenance of lymphoblast

Tissue Culture:

human lymphoblast cell line (gift of Engelward lab, MIT)

TK6 cells were grown in a flask with 12ml RPMI

supplemented with FBS. The cells were kept in an

incubator at 37°C. A stain was used to assess if the

cells were alive or dead.

10% FBS, Pep Strep X units/ml

5% CO<sub>2</sub>, 95% RH

Trypan Blue

# Improving a Methods paragraph

## Maintaining lymphoblastoid cell line(s):

- TK6 human lymphoblastoids (gift of the Engelward Lab, MIT, Cambridge MA) and Coriell cells (GM15221, GM15242, and GM15061, Coriell Institute for Medical Research, Camden NJ) were cultured at  $1-9 \times 10^5$  cells/mL in suspension in 1X RPMI medium 1640 (Invitrogen, Carlsbad CA) supplemented with 10% fetal bovine serum (Invitrogen). All cell culture media were supplemented with 100 units/mL penicillin-streptomycin (Invitrogen).

# average # cells per microwell

Yellow Team	
A1: 120,000 cells, 40 min loading	2.13, 3.02
A2: 120,000 cells, 40 min	
B1: 30,000 cells, 40 min	2.72, 2.60
B2: 30,000 cells, 20 min	
Green Team	
A1: 200,000 cells, 30min loading	2.40, 1.70
A2: 200,000 cells, 30min loading	
B1: 150,000 cells, 30min loading	1.70
B2: 200,000 cells, 15min loading	
Blue Team	
A1: 100,000 cells, 35min loading	2.3, 3.2, 1.66
A2: 100,000 cells, 35min loading	
B1: 200,000 cells, 35min loading	
B2: 100,000 cells, 15min loading	1.57
Pink Team	
A1: 40,000 cells, 30min loading	1.6
A2: 75,000 cells, 30min loading	
B1: 75,000 cells, 30min loading	1.15
B2: 75,000 cells, 5min loading	
Purple Team	
A1: 150,000 cells, 20min loading	3ish
A2: 150,000 cells, 20min loading	
B1: 15,000 cells, 20min loading	
B2: 150,000 cells, 5min loading	

loaded dead @  
~50K

lowest loading #s  
30K  
15K  
40K

# In lab today

## 1. Load TK6 to “live” chip:

- 50,000 cells/well
- For 25 min
- during cell loading incubation prepare H<sub>2</sub>O<sub>2</sub> dilutions on ice

## 2. Overlay with 1% LMP agarose

- 3 min at room temperature + 3 min at 4 °C

## 3. Treat 2 CometChips with H<sub>2</sub>O<sub>2</sub>

- one prepared by faculty, with “dead” (already lysed) cells
- one prepared by you today, with live TK6 cells

# **Thesis Defense**

Wednesday, September 21st (tomorrow!)

**Time:** 3:30 PM

**Location:** 56-614


**Investigation into the Role of DNA Damage  
and Repair during Influenza Infection and  
Inflammation**

by Marcus C. Parrish


Engelward lab

# M1D2: Analysis of Doubling Time Expt.


in main lab (T=0 days)



in TC room (~~T=4~~ 25 days)



# Doubling time of TK6 cells – sample & images



- T = 0
- fixed overnight at 4 °C with 10% formalin (on Thurs 09/15 at 5pm)
- labeled overnight at 4 °C with Vybrant Dye Cycle Green (Life Tech.)
- (not lysed)

- T = 2.5 days = 62 h
- fixed 2 h at 4 °C with 10% formalin (on Sunday 09/18 at 10am)
- labeled overnight at 4 °C with Vybrant Dye Cycle Green (Life Tech.)
- (not lysed)

VERY sensitive DNA stain



# Doubling time of TK6 cells – Matlab analysis

	A1		A2		B1		B2	
day	2.5	0	2.5	0	2.5	0	2.5	0
average signal (a.u.)	82052	4451	82052	8444	85916	10120	70377	6951
<b>doubling time (h)</b>	<b>15</b>		<b>20</b>		<b>21</b>		<b>19</b>	
<b>total DNA in microwell</b>	123495	6479	123495	1065	150278	3949	85501	9122
(a.u. = fluorescence intensity)	92015	837	92015	2843	61555	3343	85621	7244
	96638	912	96638	3802	60232	9140	43434	3925
	141691	4347	141691	1652	61232	3821	89521	8736
	141351	3949	141351	2676	48203	9981	47265	7176
	168353	3836	168353	8527	25169	2311	34547	5509
	74210	3161	74210	1834	43736	10551	32500	5140

- Fluorescence signal intensity: 1600 a.u. ~ 1 cell

$$DoublingTime = \frac{ElapsedTime}{NumberOfDoublings} = \frac{T_{end} - T_{beginning}}{\log_2\left(\frac{\#cells_{end}}{\#cells_{beginning}}\right)} = \frac{62hrs}{\log_2\left(\frac{signal_{t=2.5d}}{signal_{t=0}}\right)}$$

- Published TK6 doubling time = 20 h