M1D3: Use immunofluorescence staining to assess repair foci experiment

- 1. Prelab
- 2. Antibody staining for TDP43 localization



Overview of Mod 1 experiments

Research goal: Use functional assays to characterize ligands identified as binders to TDP43 from SMM technology



Using immunofluorescence: Localization of TDP43 in CAD cells



Using immunofluorescence: Localization of TDP43 in CAD cells

Condition 1:



TDP43

DAPI



Merge



Condition 2:



protein of interest	TDP43
primary antibody	<pre>k rabbit anti-mouse anti-TDP43</pre>
secondary antibody	人 goat anti-rabbit
Fluorophore (conjugated to secondary antibody) exc./ em. wavelengths	488/525 nm

Considerations for using antibodies in the lab

- Antibodies bind to specific epitopes on antigens
 - Antigens may have multiple epitopes
- Primary antibody recognizes the antigen
 - Specific protein sequence
 - Specific conformation of protein
 - Specific state of protein (i.e. phosphorylation)



- Secondary Ab recognizes the species of the primary Ab
 - Often conjugated to tag for visualization
 - Enzyme or fluorophore
 - Amplifies signal through multiple bindings
 - Consider sample species when choosing antibodies!



Secondary antibody conjugated to a fluorophore

Epițopes

Polyclonal vs. monoclonal antibodies

Polyclonal

- How it's made: animal (often rabbit) immunized with antigen of interest then antibodies collected from blood sera and affinity purified
- Advantages:
 - Less expensive and faster to produce than monoclonal
 - Multiple antibodies in one polyclonal mixture can increase antigen recognition by binding multiple epitopes
 - Especially useful for proteins with low expression
- Disadvantages:
 - Variability from lot to lot



Monoclonal

- How it's made: animal (usually mouse) immunized with antigen of interest then B cells from spleen are harvested and fused with myeloma cells to create hybridoma cell line that will continually produce single antibody clone
- Advantages:
 - Very consistent
 - Binds single epitope (can also be disadvantage)
- Disadvantages:
 - More expensive and requires animal sacrifice



Using immunofluorescence (IF): steps in protocol

ParaFormaldehyde Fix cells



Finish IF by adding DAPI, then mount slides for imaging



Mount coverslip on glass slide with mounting media

Blue= DAPI Green= antibody staining

For today:

- 1. Complete IF staining for TDP43 Localization
 - 1. Christine will demo staining chamber assembly
- 2. Work on Methods revision with partner

For M1D8

- 1. <u>Individually</u>, answer the question prompts for the Implications and Future works section of your Data Summary
- 2. <u>With your lab partner</u>, revise your methods homework and add M1D4-M1D5