

M1D2: Complete small molecule microarray analysis and induce protein expression

Announcements

Office hours:

Mon 2-5pm (Noreen, 16-317)


Tu 4-5pm (Josephine, 56-341c/322)

Th 10-11am (Josephine, 56-341c)

Email us for other times

1. Pre-lab discussion
2. Induce FKBP12 expression
3. Gel electrophoresis confirmation digests
4. Complete SMM data analysis

Overview of Mod1 experiments




Review of cloning: generate plasmid to make FKBP12 protein


1) Amplification

2) Digestion

3) Ligation



A closer look at the pRSETb expression vector (a.k.a. backbone)



- P_{T7} promoter
engineered to be IPTG-inducible
- 6xHis
6 histidine tag
binds to Nickel
Helps with purification
- MCS
Multiple cloning site

M1D1 Part 2c: Calculate volumes of insert and backbone needed for ligation

Knowns:

- Need 50-100 ng backbone (specified by manufacturer)
- pRSET vector: 2853 bp
- FKBP12 Insert: 331 bp
- Molar mass $\sim 660\text{g}/(\text{mol} \cdot \text{bp})$
- Desired molar ratio of insert to backbone is 4:1
- Concentration of backbone (from recovery gel):
- Concentration of insert (from recovery gel):

Calculate 4:1 (insert:backbone) *molar* amounts, final volumes for ligation

1. Calculate moles of backbone (50-100 ng, e.g. choose 60 ng)

- $2853 \text{ bp} * (660 \text{ g} / (\text{mol} * \text{bp})) = 1.88 \times 10^6 \text{ g/mol}$
- so $60 \text{ ng} / (1.88 \times 10^6 \text{ g/mol}) = 3.19 \times 10^{-14} \text{ mol}$

2. Determine moles of insert needed (4x backbone)

- $4 \times 3.19 \times 10^{-14} \text{ mol} \sim 1.26 \times 10^{-13} \text{ mol}$
- with $331 \text{ bp} * (660 \text{ g} / (\text{mol} * \text{bp})) = 2.18 \times 10^5 \text{ g/mol}$
- so use $1.26 \times 10^{-13} \text{ mol} * 2.18 \times 10^5 \text{ g/mol} \sim 27 \text{ ng}$

3. Calculate volume of backbone and insert needed

- Backbone: $60 \text{ ng} / (20 \text{ ng}/\mu\text{L}) = 3 \mu\text{L}$
- Insert: $27 \text{ ng} / (40 \text{ ng}/\mu\text{L}) = 0.68 \mu\text{L}$

MIDI Part 2C: Ligation Calculation

Recovery gel:

backbone: $100\text{ng}/5\mu\text{L} = 20\text{ng}/\mu\text{L}$

insert: $200\text{ng}/5\mu\text{L} = 40\text{ng}/\mu\text{L}$

① Volume of backbone (50-100ng)
eg. let's choose 60ng \rightarrow 3 μL

② Calculate moles of backbone

$$\text{backbone} \quad \frac{60 \times 10^{-9} \text{ g}}{660 \text{ g}} \cdot \frac{\text{mol} \cdot \text{bp}}{2853 \text{ bp}} =$$

$$3.19 \times 10^{-14} \text{ moles}$$

③ Determine moles of insert needed
(4:1 molar ratio of Insert to backbone)


$$4 \times 3.19 \times 10^{-14} \text{ moles} = 1.26 \times 10^{-13} \text{ moles}$$

④ Calculate volume of insert needed

$$1.26 \times 10^{-13} \text{ mol insert} \left(\frac{660 \text{ g}}{\text{mol} \cdot \text{bp}} \right) \left(\frac{3316 \text{ bp}}{1} \right) \left(\frac{\mu\text{L}}{4 \times 10^{-8} \text{ g}} \right)$$
$$= 0.68 \mu\text{L insert}$$

How do we confirm the plasmid product?

- Amplify plasmid
 - Transform into bacteria
- Purification
 - Separate plasmid from chromosomal DNA
- Digestion [This is what you did last time \(M1D1\)](#)
 - Confirm the plasmid contains expected fragments




Gel electrophoresis to visualize confirmation digest results

Driving force:
electric field,
charge

DNA is
+ / - charged

Separates DNA
by size



1% agarose gel

Visualize DNA + save a picture!

- DNA Loading dye (6X):
 - bromophenol blue**
 - runs as ~500 bp**
 - visualize progress**

glycerol

-viscous, dense

-DNA sinks in well



- Sybr-Safe DNA stain:

-DNA intercalator


-Fluorescent

-Visualize with UV or blue light


- **Safety : wear nitrile gloves**



Your colleagues in Spring 2017...



Using immunofluorescence to detect ligand binding





SMM quantification steps

1. Align GAL file to fluorescence on 532 nm channel (sentinel spots)
2. Quantify fluorescence on 635 nm channel
3. Identify 'hits' with improbably high fluorescence
4. Identify compounds that hit repeatedly
5. Compare top hits to common binders list

Align SMM results to using sentinel spots

- Every spot can be located using intersecting lines between sentinels




Images represent arrays of numbers

- Each pixel is a 16-bit number that represents intensity
- Computational analysis used to define 'hits'

4	3	4	4	3	2	3	4	3	5	4	6	3	3	3	2	3	2	2
3	5	4	3	3	3	5	6	7	8	5	6	4	4	4	3	3	3	3
3	3	3	3	4	8	12	92	275	311	256	61	11	6	3	3	3	3	4
4	3	3	4	8	173	625	818	823	856	815	831	568	136	9	5	4	4	3
5	3	4	8	273	830	814	835	873	890	836	857	818	771	201	9	6	2	2
3	4	7	175	780	805	877	941	936	920	973	921	842	819	714	125	6	3	2
4	4	29	568	868	867	905	909	936	994	954	931	963	875	813	490	15	5	4
4	5	131	754	852	906	958	920	963	923	917	904	951	930	851	716	95	6	3
4	5	229	796	879	924	934	923	962	961	993	993	945	989	867	780	162	6	4
3	7	254	827	879	965	949	960	982	926	918	955	927	984	872	765	204	7	3
4	5	175	808	883	996	951	998	935	976	971	940	922	961	872	804	132	4	4
4	4	57	666	859	968	999	947	977	985	916	928	960	974	841	678	62	4	4
4	3	11	406	839	897	915	930	946	993	914	911	977	900	830	359	10	3	4
3	2	5	60	624	830	890	973	903	921	912	930	881	850	613	54	6	3	3
3	4	4	7	92	602	873	856	882	913	887	885	842	589	82	7	4	3	3
3	4	3	4	5	23	266	697	838	828	837	667	261	21	5	4	4	5	4
3	3	4	4	4	6	9	12	27	49	28	11	9	7	5	3	3	4	3
3	5	3	5	4	4	7	4	4	6	6	3	5	3	3	3	3	4	4

Quantify fluorescence to identify hits

- Foreground
- Background



$$\text{Signal-to-noise ratio (SNR)} = \frac{\mu_{\text{foreground}} - \mu_{\text{background}}}{\sigma_{\text{background}}}$$

Identify hits based on z-score (typically is a measure of how many standard deviations away from the mean the signal is). We will use a robust z-score, which is the same idea but based on the median intensity instead of the mean.


signal

noise

How to evaluate the SMM results

- Is the background noisy?
- Are the positive controls easily recognized?
- Do any areas appear strange? Damaged?
 - Manufacturer or handling defects
- Are the hits aligned with printing spots?
- Do you trust the data?

What do we expect to see?



Factors that influence hit identification

- How many false positives are expected?
 - More hits needed if confidence is low
- How many chemical ‘patterns’ are evident?
 - Repeated patterns between compounds may increase confidence
- Are the hits unique to the screen?
 - Promiscuous binders may decrease confidence

Today in lab...

- Reminder: wipe down bench with 70% EtOH before and after wetlab work
- Reminder: empty benchtop waste bucket into biowaste bin at end of day
- Measure OD of bacterial culture—when OD is >0.5 , add IPTG to induce expression of FKBP12 protein
- Perform electrophoresis on confirmation digest reactions
- SMM data analysis

For next time...

- Draft a figure with your confirmation digest results for your Data Summary
 - Include a title and caption
- Schedule appointment at BE Communication Lab before M1D5

Notes on figure making:

- Image should not be the entire page
 - Only needs to be large enough to be clear
- Title should be conclusive
 - Don't include what you did, rather include what you found
- Caption should not detail the methods
 - Define abbreviations, symbols, etc.