


## M1D2: Prepare and treat cells for $\gamma$ -H2AX experiment

1. Prelab

2. Treat and fix MEF cells for H2AX assay

3. Work on Background and Motivation section



## What are we studying in Mod 1 and why?

Research question: Does exposure to As inhibit repair of H<sub>2</sub>O<sub>2</sub>induced DNA damage, raising the possibility that combined exposure
is an important risk to public health?

Why are we interested in DNA Damage?

Why are we interested in arsenic?

• Why would we use  $H_2O_2$  to induce DNA Damage?

## Mod1 Overview (AKA: how we approach our research question)

Last lab: Seeded

MEF cells on

Coverships For


YHZAY assay


This lab:

Treat alls w/ As & H202 Fix cells

Next lab:

coverslips >>
Ab staming
YIFZAX





- 1. Use repair foci experiment to measure DNA breaks
- Examine effect of  $H_2O_2$  +/- As on double strand DNA breaks by measuring  $\gamma H2AX$  foci formation

DSB-double Strand whats

- 2. Use high-throughput genome damage assay to measure DNA damage
- Measure effects of H<sub>2</sub>O<sub>2</sub> +/- As on DNA damage by measuring DNA migration in agarose matrix

SSB = Single Strand broaks

# How does $H_2O_2$ damage DNA? H<sub>2</sub>O<sub>2</sub> + e<sup>-</sup> $\rightarrow$ HO<sup>-</sup> + •OH — hydroxyl

$$- \rightarrow HO^- + \bullet OH - nydnxyl$$

$$2 H^+ + 2 e^- + H_2O_2 \rightarrow 2 H_2O$$

By producing Reactive Oxygen Species

8-oxo-Guanine


## How does H<sub>2</sub>O<sub>2</sub> damage DNA?

#### **Mutation if replicated**



## How do cells repair oxidative DNA Damage?

Typically, the BER pathway



## How do we look at DNA damage in intact cells?


#### Look for γ -H2AX foci

H<sub>2</sub>O<sub>2</sub> can cause damage to DNA, resulting in a damaged base

That damage causes a single strand break as cell tries to repair the DNA

Multiple single strand breaks cause double strand breaks


At DNA strand breaks, DNA-PK and ATM phosphorylate the histone H2AX




We can identify the frequency of these breaks as a measure of DNA damage

## Poly-d-Lysine

- Want clear images  $\gamma$  -H2AX foci in nuclei
  - Immobilize cells in a monolayer on glass coverslips
  - Even adherent cells can benefit from a coated surface



- Many ways to get cells to adhere to glass or plastic
  - ECM molecules (like laminin)
  - Charged polymer molecules (like PdL)
  - Biological substrates (like Polyphenolic Proteins secreted by marine mussels)



## Treatment conditions for this experiment

- Goal: identify any additive effect pretreatment with As has on H<sub>2</sub>O<sub>2</sub> induced DNA damage
  - Treat cells with As for 4 hours, then treat cells with H<sub>2</sub>O<sub>2</sub> for 15 minutes


**Experimental Condition** 

tooz + AS Dum 10em 40mm **Control Conditions** 

As above H202 alone no treatment

## Fixing cells with paraformaldehyde

- Formaldehyde is a chemical crosslinker
  - Paraformaldehyde is the polymer version depolymerize into formaldehyde with heat or basic conditions
- Commonly used in cell and tissue microscopy
  - Crosslinking of proteins will "fix" samples in place for static imaging
  - Can also trap nucleic acids, etc... between crosslinked proteins



Kiernan, 2000

## Homework discussion

Methods

## Pro tips for writing a methods section

#### Include enough information to replicate the experiment

- List manufacturer's name (Company)
- Be concise and clear in your description

#### Use subsections with descriptive titles

Put in logical order, rather than chronological order
 Begin with topic sentence to introduce purpose / goal of each experimental procedure

#### Use clear and concise full sentences

- NO tables or lists, all information should be provided in full sentences and paragraphs
   Write in passive voice and use past tense

#### Use the most flexible units

Write concentrations (when known) rather than volumes

#### Eliminate 20.109 specific details

- Example "labeled Row A, Row B..."
- Do not include details about tubes and water!
- Assume reader has some biology experience

Include parts of the protocol that the teaching faculty completed, but do not say "completed by teaching faculty."





## How can you improve this example?

"Cells were grown in 12 mL of DMEM supplemented with FBS. We spun

down the cells and counted them with a hemocytometer.

Is know long, speed, other conditions

were incubated in 37 C incubator."

### How can you improve this example?

What cells? From where were the cells attained?

How much? What else was added to the media?

"Cells were grown in 12 mL of DMEM supplemented with FBS. We spun

Volume here does not have context as based on the flask used. When might flask / plate size be helpful?? Define all abbreviations and include supplier / manufacturer.

Use passive voice and avoid jargon!

down the cells and counted them with a hemocytometer. Flasks

Be specific about the purpose of each of the steps used...cells were harvested using centrifugation (be sure to include speed and time) then counted using a hemocytometer. And what else was used? At what final concentration / percent?

Be specific about the subject of each action / step.

were incubated in 37 C incubator."

Specific location / equipment used is not important, just the temperature conditions. What other growth conditions were maintained?

## Revised example...

J.St.

Maintaining MEF cell line

Mouse Embryonic Fibroblasts (MEF) cells (gift of Engelward Laboratory, MIT) were grown in Dulbecco's... (DMEM) (Manufacturer) supplemented with 10% fetal bovine serum (FBS) (Manufacturer) and 100 U / mL of penicillin and streptomycin (Manufacturer). To harvest, cells were trypsinized, collected in fresh media, and counted using 10% (v/v) trypan blue and a hemocytometer. Cultures were maintained at 37 C, 5% CO<sub>2</sub>, and 95% relative humidity.

## In lab today

- 1. Expose As-treated cells to H<sub>2</sub>O<sub>2</sub> and fix them for IF staining
- Talk about purpose and structure of background and motivation following the experiment

DO NOT use the aspirator to remove cell media with Arsenic. It needs to be collected as hazardous waste!

#### M1D3HW

- 1. Work with your lab partner to write methods for M1D1 and M1D2
- 2. Schedule appointment with BE Comm Fellow before M1D5

What is the purpose of the Background & Motivation section?

What *specific* information should be included in the Background & Motivation section?

## Notes on Background & Motivation section...

- Anchor your research in a general topic that is important to a broad audience
  - Focus on describing what is currently known in the field
  - Reference the relevant research in the field
- Connect your research to the general topic
  - Minimum essential information
  - Introduce specific technologies necessary for understanding your specific project
- Address how you will expand on what is currently known
  - Include evidence of incompleteness of current understanding
  - Motivate your investigation
  - Include a clear hypothesis / research goal
- Provide a preview of your findings and the implications
  - Tie back to the initial general topic
  - Avoid including extensive methods details

## Notes on topic sentences...

 Topic sentence = first sentence of each paragraph

- Should 'funnel' from big picture topic to your specific research question / project
  - Provide only the background needed to understand research / problem / goal
  - Clearly state what is not currently known
  - Address how you will fill knowledge gap
  - Provide preview of your results

• Include references!!

**Impact Statement** 

Specific background

Knowledge gap/ Statement of problem

Hypothesis

Here we show...

## How should you introduce your story?

1<sup>st</sup> paragraph: what is the big picture / problem?

2<sup>nd</sup> paragraph: what is currently known?

3<sup>rd</sup> (or 4<sup>th</sup>) paragraph: what is your research question?

4<sup>th</sup> (or 3<sup>rd</sup>) paragraph: how will you address your question?

5<sup>th</sup> paragraph: here we show...