

# M2D3: Purify RNA and practice

## RNA-seq data analysis methods

03/15/2018

1. Prelab discussion
2.  $\frac{1}{2}$  class to TC to harvest cells for RNA purification followed by cDNA synthesis
3.  $\frac{1}{2}$  group start practice RNA-seq analysis

# Mod2 major assignments


- **Research Article (20%)**
  - individual, submit on Stellar
  - due April 21<sup>st</sup> at 10pm
  - format: word document
- **Journal Club Presentation (15%)**
  - individual, presentation during lab
  - presentation slides due on Stellar 1pm April 3<sup>rd</sup> or April 5<sup>th</sup>
  - format: powerpoint, keynote, or google slides
- *Lab quizzes (5%)*
- *Homework and Notebook (10%)*
- *Blog (5%)*
  - *by Sunday, March 18 at 10 pm (Mod1)*
  - *by Saturday, April 7 at 10 pm*
  - *by Sunday, April 22 at 10 pm*

# 20.109(S18) Class blog


- Possible topics listed on the blog website
- Details about use:
  - Do not publish MIT logo
  - Do not post photographs with names tagged
  - Do not write malicious comments
  - Do not plagiarize



# M2: Experimental overview



**Analyze RNA-seq results**




DAY 4: Evaluate altered gene expression


DAY 5: Investigate public databases

DAY 8: Identify regulatory motifs

**Examine effect of drug treatment on cell survival**



# M2D2: (Noreen)Treated cells with etoposide



M2D3: extract RNA (~48 hours after DNA damage)

# Qiagen Isolate RNA: QIAshredder + Rneasy kit



purple



pink



steps	contents	purpose
lyse	RLT (with highly denaturing guanidine-thiocyanate salt) + QIAshredder	inactivate RNase, disrupt membranes, helps bind column  homogenize (shear high-MW genomic DNA)
prepare	70% ethanol	promote efficient binding to silica
bind	silica membrane in column	retain mRNA
wash	RW1 RPE	remove contaminates ** after this wash, important to get rid of <u>all</u> ethanol
elute	water, RNase-free	high-purity RNA

**RLT buffer:** composed of detergents and chaotropic salts(weakens hydrophobic effects)

**Qiasredder:** polymer that shears high molecular weight components of the cell

**EtOH:** RNA insoluble in ethanol, RNA precipitates from cell lysate and binds to silica membrane

# Components and procedure of cDNA Synthesis

steps	conditions	reagents added
denature & anneal	65°C 5 min on ice 1 min	<u>1 µg RNA</u> + oligo (dT) <sub>20</sub> primer + dNTPs (dATP, dCTP, dGTP and dTTP)
synthesize cDNA	50°C 50 min	<u>Superscript III Reverse Transcriptase</u> MgCl <sub>2</sub> DTT RNase OUT buffer
terminate	85°C 5 min	<i>kills RT enzyme</i>
remove RNA	37°C 20 min	RNase H
Purify cDNA	M2D5	

What genes are differentially expressed in response to DNA damage?

*How are we addressing this question?*

DLD-1

DLD-1  
+  
etoposide

BRCA2-/-

BRCA2-/-  
+  
etoposide

purify RNA

↓  
synthesize cDNA


↓  
quantitative PCR: PC1, GAPDH  
(2 genes)

→ Illumina  
RNA  
sequencing  
↓  
all genes




# Workflow for Illumina HiSeq 2000


(1) sample prep



(2) amplification



(3) sequencing



(4) data analysis



# Reminders:

- M2D4 HW: Methods M2D1-M2D3 (omit practice RNA-seq analysis)
  - NO use of “per manufacturers protocol” for this methods homework
- Mini presentation due Saturday March 17<sup>th</sup> at 10pm.
  - Email video file to [bioeng20.109@gmail.com](mailto:bioeng20.109@gmail.com)
  - Submitting the final version of your video can take time so don't wait till the last minute. Feel free to send us a link so we can download.