


M1D2:


Perform protein purification protocol

1. Prelab discussion
2. Purify MAX-6xHis protein
3. Electrophoresis confirmation digest



“Don’t pick it up,” I say, and he says, “It’s just a *plasmid*, what harm could it do?” Well just look at him now...who knows *what* protein he’s expressing!

Overview of Mod 1 experiments:




What is our protein of interest?

- MAX functions as a transcription factor
 - Forms homodimers and heterodimers
- Dimerizes with Myc, which is an oncogenic transcription factor
- Homodimers and heterodimers compete for binding at promoters to provide regulatory system of target genes



Closer look at pET28a_MAX-6xHis



The MAX-6xHis insert was cloned into the pET28a expression vector. This vector was engineered for over-expression of proteins via the incorporation of the components listed on the right.

- 6xHis
- *lac* operator
- *lacI*
- T7 promoter
- T7 terminator
- RBS

'Tag' that will allow us to purify the expressed MAX protein

Sequence found in bacterial genomes that is used to control transcription of certain genes


Gene that encodes LacI, repressor protein that inhibits transcription of certain genes

Promoter found in phage that is used to control certain genes

Terminator found in phage that is used stop transcription

Ribosome binding site sequence


Overview of protein expression system




BL21 is a specialized strain of E.coli designed to make protein

- 1) T7 RNA Polymerase is added
- 2) Protease genes removed

T7 RNA polymerase transcribes MAX-6xHis




T7 RNA polymerase transcribes MAX-6xHis



The *lac* operon evolved to induce context-sensitive transcription of lactose metabolism genes – we add a lactose biomimic!

How did we induce protein expression?




These steps were completed by the instructors for timing reasons. Please note that this information should still be included in your methods homework!


How will you purify MAX-6xHis?

- First, need to lyse cells to release proteins
 - B-PER (Bacterial Protein Extraction Reagent):
Detergent in sodium phosphate buffer. Aids in lysis by solubilizing lipids and proteins in the membrane, thereby creating pores.
 - Lysozyme:
Enzyme solution that contains lysozyme and benzonase. Lysozyme cleaves glycosidic bonds in peptidoglycan and benzonase cleaves nucleic acids.
 - Proteinase inhibitor:
Molecules that bind the proteolytic (active) sites of proteinase enzymes.
 - Dnase:
Cleaves DNA. Added in addition to the benzonase as pilot experiments demonstrated that the DNA was not completely cleaved. Uncleaved DNA is quite viscous and can clog the column during protein purification.


Overview of protein purification steps used after lysis of cells:



6xHis tag binds to Ni²⁺ resin / column





- Ni²⁺ chelated onto agarose resin via nitrilotriacetic acid (NTA) ligand



- His tag chelates to Ni²⁺ causing protein to 'stick' to resin / column

The basic imidazole ring of histidine binds to the Ni²⁺ ions causing tagged MAX protein to bind to the agarose resin while other proteins pass through the column.


Non-specific binders washed from Ni²⁺ resin / column using imidazole



- Low concentration of imidazole included in wash buffer
- Imidazole competes for binding to Ni²⁺ resin
 - Low affinity binders / non-specific binders are outcompeted and released from the resin


Because other proteins in the cell have histidine residues, it is possible to have non-specific binding to the column. A low concentration imidazole solution is used to wash weakly bound proteins from the column.

Imidazole used to elute protein from column




- Elution buffer contains higher concentration of imidazole compared to wash buffer

Binding:




Elution:



- Increased concentration allows imidazole to out-compete 6xHis for binding to Ni^{2+} resin

This takes advantage of the chemistry that enables the His tag to bind to the Ni^{2+} ions. The imidazole ring of His binds MAX-6xHis to the column. This interaction can be out-competed by adding a high concentration imidazole solution to the column.

Wrap-up of confirmation digest



- How do you visualize migration through the gel?

Loading dye is used to visualize migration during electrophoresis. The components are bromophenol blue (dye) and glycerol. The glycerol 'weights' the sample down into the well so it doesn't 'float' into the buffer.

- How do you visualize DNA bands in the gel?

SYBER safe is used to visualize the DNA bands after electrophoresis. This is a fluorescent dye that binds to DNA. The DNA-dye complex absorbs blue light and emits green light.

For today...

- Watch Protein purification demo
- Start protein purification protocol
- Complete agarose gel electrophoresis during lysis incubation
- **Be sure to clearly label all tubes containing protein purification aliquots!**

For M1D3...

- Draft methods section for confirmation digest (M1D1) and protein purification (M1D2) protocols.

Notes on methods section...

Include **enough information to replicate** the experiment

- Cite manufacturer for supplies / equipment (Company)
- Be concise and clear in your description

Use subsections with **descriptive titles**

- Put in logical order, rather than chronological order
- Begin with topic sentence to introduce purpose / goal of each experimental procedure

Use clear and concise **full sentences**

- NO tables or lists, all information should be provided in full sentences and paragraphs
- Write in passive voice and use past tense

Use the most **flexible units**

- Write concentrations (when known) rather than volumes

Eliminate 20.109 specific language and obvious details

- Example “labeled Row A, Row B...”
- Do not include details about tubes and water!
- Assume reader has some biology experience
- Include parts of the protocol that the teaching faculty completed, but do not say “completed by teaching faculty.”

How can you improve this example?

DNA was cut to check insert. Enzymes were used for single and double

digest then run on gel made by adding 1 g of agar to 100 mL of water.

Gel was imaged on a gel box.

How can you improve this example?

What DNA?
How much?

Consider more formal language.

What insert?

Which enzymes? From where were the enzymes acquired?

DNA was cut to check insert. Enzymes were used to cut DNA for

Specifically, why was this done?

Redundant

Provide details on how this was done.

Colloquial...use more scientific language. Also, include details.

single and double digests then run on a gel made by adding

What does this mean?

Be mindful of the order of information and of confusing sentence structure.

1 g of agar to 100 mL of water. Gel was imaged on the gel box.

Use the most flexible units / concise description.

What else was needed for imaging?

The?

What would be more informative?

Edited example...

Confirmation digest of pET28a_MAX-6xHis

To confirm that MAX-6xHis was cloned into pET28a expression vector, a digest was completed. Restriction enzymes Abcl and Defl were used to digest W ng of pET28a-MAX-6xHis in single digests (only one enzyme added) and in a double digest (both enzymes added) using X U / uL of each enzyme and Y buffer (NEB). Digests were incubated at 37C for Z hrs then stored at -20C. [Add gel electrophoresis details].

You can use this as part of your homework assignment! Be sure to include the details for placeholders and the information for the gel electrophoresis protocol.