

M1D6: Examine sub-nuclear foci abundance to measure DNA damage

10/03/17

1. Pre-lab Discussion
2. Start H2AX staining of fixed cells
3. Process CometChip microscopy images and analyze via Matlab
4. Rotate groups to Engelward microscope to image CometChip
5. 4:15pm start paper discussion

Assignments for M1

- Data summary draft
 - due by 10pm on Wed., October 11
 - revision due by 10pm on Sun., October 22

Summary content

1. Title
2. Abstract
3. Background, Motivation
4. Figures, Results & Discussion, Interpretation
5. Implications, Future Work

- Mini presentation due by 10pm on Sat., October 14
- Blog post for M1 due by 10pm on Mon., October 23

Extra office hours


- **56-302**
- Saturday, 10/7, 10am-2pm
- Wednesday, 10/9, 10am-3pm

Regular office hours

- **Next week Monday OH canceled**
- Josephine, Thursday 2pm-3pm in 56-341c
- Leslie, Friday 9am-10am and 3pm-4pm in 56-341c


Overview of

Module 1: Measuring Genomic Instability




1. Optimize comet chip assay

- Test loading variables



2. Use comet chip assay to measure DNA repair

- Measure effects of MMS and H_2O_2 on BER
- Measure (Quantify) effect of DNA repair activity in response to alkylation damage and oxidative damage.





3. Use immuno-fluorescence assay to measure DNA repair

- Examine effect of MMS and H_2O_2 on DSB abundance
- Quantify effect of DNA damage response (DNA repair and recovery) in response to alkylation damage and oxidative damage induced DSBs.

recovery= continued repair w/o presence of DNA damaging agents


enzyme treatment minus the non enzyme treatment= actual (unrepaired)
damage from MMS/H₂O₂

What result do we expect from our
assay so far...



Incubating the Comet Chip with purified enzymes
(OGG1, AAG) **REVEALS** damage not quantified
otherwise (buffer control chip)


Base Excision Repair (BER) Pathway



not necessarily visualized in comet tail

strand breaks visualized via Comet tail



Complete experiment to test role of biochemical factors (MMS or H2O2) in genomic stability (DNA damage)



Measuring Double strand breaks via γ H2AX Assay




Fix cells and stain with
antibody that marks
double strand breaks








Blue: DNA
Green: γ H2AX staining=
double strand breaks

Using immunofluorescence: γ H2AX assay to detect double-strand DNA breaks




- histone H2AX phosphorylated at Ser139 if DSB
- use antibodies against γ H2AX




protein of interest	 γ H2AX
primary antibody	 mouse anti-human anti- γ H2AX
secondary antibody	 goat anti-mouse
fluorescent dye exc./ em. wavelengths	 488 / 520 nm

M1D6, practically using immunofluorescence: γ H2AX assay to detect double-strand DNA breaks



What do we hypothesize will happen...



Cells fixed @30min
No recovery #1

Cells treated for 30min
Recovery in fresh media #2

CometChip Data analysis in ImageJ and MATLAB

1. ImageJ


- from several images per well to one stack per well
- GenImageStacks_singleimage.txt

2. MATLAB


- for each comet in stack, calculates intensity of head and tail, as well as length of tail

3. Excel

- export data from Matlab and compile
- post data to the wiki



What's in the final Excel file?



	01	02	03	04	05	06
A	7.45	7.68	11.33	16.49	34.06	29.43
B	8.59	7.33	10.03	14.49	26.58	37.04
C	6.86	8.73	11.94	18.78	34.69	37.87
D	10.37	11.93	10.77	12.14	9.68	11.71
E	14.10	10.54	9.76	10.79	11.85	10.32
F	15.28	10.51	9.53	10.36	11.67	9.29
[H2O2] (mM)	0	0.25	0.5	1	2	4
[MMS] (uM)	0	10	20	40	60	80

triplicates

Cometnumbers	%Head DNA	%Tail DNA	OTM (um)	Tail Len. (um)	Comet Len. (um)	+
--------------	-----------	-----------	----------	----------------	-----------------	---

- Cometnumbers: how many comets were used for calculation in each well (= stack)
- %Head DNA = $100 * \text{HeadFluorescence} / (\text{HeadFluorescence} + \text{TailFluorescence})$
- %TailDNA = $100 * \text{TailFluorescence} / (\text{HeadFluorescence} + \text{TailFluorescence})$
- Olive tail moment (OTM) = $(\% \text{TailDNA} / 100) * (\text{TailCenterOfMass} - \text{HeadCenterOfMass})$
- Tail length
- Comet length

Don't forget reported data must have 95% C.I. and include statistical analysis

- Plot your data with 95% confidence intervals

$$\bar{x} \pm \frac{t_{table} * stdev}{\sqrt{n}}$$

$$t_{table} = TINV(0.05, n - 1)$$

- How certain are you that two populations are different?

$$p = TTEST(array1, array2, 2, 3)$$

2-tailed

unequal variance (heteroscedastic ☺)

- ✧ The Student's t-test only applies to **two** data sets.

Only compare two conditions at a time.

In lab today

1. Retrieve your 12 well plates and choose 6 conditions to for your immunofluorescence experiment; start staining protocol.
2. During incubations begin Matlab analysis.
3. Groups will go one at a time to the Engelward microscope to see your CometChip.
4. Paper discussion will start at 4:15pm in the lab.

Cytometry

PART A
Journal of the
International Society for
Advancement of Cytometry

Standard Fluorescent Imaging of Live Cells is Highly Genotoxic

Jing Ge,¹ David K. Wood,² David M. Weingeist,¹ Somsak Prasongtanakij,³ Panida Navasumrit,³ Mathuros Ruchirawat,³ Bevin P. Engelward^{1*}

HWM1D7 Methods M1D1-M1D5

- Complete assignment as team
- First sentence in each subsection should be a brief introductory sentence motivating the method
- Use final concentrations
- Don't '109 specific language
- Be concise while still providing enough information for the reader to repeat the experiment
 - Example: buffer recipes and wash times/number in parentheses