



M2D2: Design gRNA sequence for CRISPRi system



- Prelab discussion
- Diagnostic digest electrophoresis
- Design gRNA primer to target (and inhibit) gene in the metabolic pathway of E. Coli



Quick Recap



CRISPRi



Mod2 Overview

Research goal: Increase the yield of commercially valuable byproducts in *E.coli* using CRISPRi technology to target genes involved in mixed-acid fermentation pathway.

Last Lab:


Cloned pdCas9 construct to generate 1st component of CRISPRi system

This Lab:


Choosing metabolic gene of interest and designing sgRNA to target it

Next Lab:

Will clone sgRNA into vector to create 2nd component of CRISPRi system



Using CRISPRi as an experimental tool



- sgRNA (**single guide RNA**) binds to region of gene of interest
 - contains a dCas9 handle
- dCas9 uses sgRNA as guide and sits on DNA without cleaving it
 - blocks transcription

Today's goal:


Design gRNA sequence to repress a gene, such that the production of ethanol or acetate will increase.

Considerations for successful gRNA primer design


Primer size range: 20-25bps

- 1) Pick a gene target to inhibit
- 2) Where is the gene?
- 3) Determine target sgRNA sequence location
 - 5' UTR?
 - Upstream regulatory element?
 - Coding region?
- 4) Template vs Nontemplate targeting
- 5) Check RNA specificity (BLAST)
- 6) Assemble sgRNA sequence

Which part of your selected gene is best to target?






Block transcription initiation





Effective for both NT and T strands

Block transcription elongation




Effective only for the NT strand

Which part of your selected gene is best to target?






Design of the sgRNA



Design of the sgRNA

Protospacer Adjacent Motif




Design of gRNA for CRISPRi system

(1) Target the **TEMPLATE** DNA strand:

gRNA sequence will be the same as the transcribed (nontemplate) sequence.

(2) Target the **NONTEMPLATE** strand:

gRNA sequence will be the reverse-complement of the transcribed sequence.




Design of gRNA for CRISPRi system

(1) Target the **TEMPLATE** DNA strand:

gRNA sequence will be the same as the transcribed (nontemplate) sequence.

(2) Target the **NONTEMPLATE** strand:

gRNA sequence will be the reverse-complement of the transcribed sequence.




Use the KEGG database to find sequence information for your target gene

Position	1937649..1939091 Genome map
AA seq	480 aa AA seq DB search MSRRLRRTKIVTTLGSPATDRDNNLEKVIAAGANVVRMNF SHGSPEDHKMRADKVREIAAK LGRHVAILGDLQGPKIRVSTFKEGKVFLNIGDKFLLDANLGKGEKDKEKVGIDYKGLPAD VVPGDILLDDGRVQLKVLEVQGMKVTFTEVTVGGPLSNNKGINKLGGGLSAEALTEKDKA DIKTAALIGVDYLAVSFPRCGEDLNYARRLARDAGCDAKIVAKVERAEAVCSQDAMDDII LASDVVMVARGDLGVEIGDPELVGIQKALIRRARQLNRAVITATQMMESMITNPMPTRAE VMDVANAVLDGTDVAVMLSAETAAGQYPSETVAAMARVCLGAEKIP SINVSKHRLDVQFDN VEEAIAMSAMYAANHLKGVTAIITMTESGRTALMTSRISGLPIFAMSRHERTLNLTALY RGVTPVHFDSANDGVAAASEAVNLLRDKGYLMSGDLVIVTQGDVMSTVGSTNTTRILTVE
NT seq	1443 nt NT seq +upstream <input type="text" value="0"/> nt +downstream <input type="text" value="0"/> nt atgtccagaaggcttcgcagacaacaaaatcggtaccacgttcagccagcaacagatcgc gataataatcttgaaaaagttatcgcggcggttgccaacgttgtagctatgcttttct cacggctcgcctgaagatcacaaaatgcgcgcggataaagttcgtgagattgccgcaaa ctggggcgctcatgtggctattctgggtgacctccagggggccaaaatccgtgtatccacc tttaaagaaggcaaggttttccatattggggataaattcctgctcgacgccaaacctg ggtaaagggtgaaggcgacaaagaaaaagtcggtatcgactacaaaggcctgctgctgac gtcgtgacctggtgacatcctgctgctggacgatggctcgcgtccagttaaaagtactggaa gttcagggcatgaaagtgttcaccgaagtcaccgtcgggtgggtcccctctccaacaataaa ggatcaacaaacttgccggcggtttgtcggctgaagcgtgaccgaaaaagacaaagca gacattaagactgcggcggttgattggcgtagattacctggctgtctccttcccacgctgt ggcgaagatctgaactatgcccgtcgcctggcacgcgatgcaggatgtgatgcgaaaatt gttgccaaggttgaaacgtgcggaagccgtttgcagccaggatgcaatggatgacatcatc ctcgccctctgacgtggtaatggttgacgtggcgacctcgggtgtggaaattggcgacctg gaactggctcggcattcagaaagcgttgatccgtcgtgcgcgtcagctaaaccgagcggt atcacggcgacctcagatgatggagtcgaatgattactaaccgatgccgacgcgtgcagaa gtcatggacgtagcaaacgcgttctggatggtactgacgctgtgatgctgtctgcagaa actgccgctgggcagtatccgtcagaaaccgttgacgcatggcgcgcggttgccctgggt gcggaaaaaatcccagcatcaacgtttcttaaacaccgtctggacgttcagttcgacaat gtggaagaagctattgccatgtcagcaatgtacgcagctaacccacctgaaaggcggtacg gcgatcatcaccatgaccgaatcgggtcgtaccgcgctgatgacctcccgtatcagctct ggctctgccaatcttcgccatgtcgcgccatgaacgtacgctgaacctgactgctctctat cgtggcggttacgccggtgcactttgatagcgctaataacggcgtagcagctgccagcgaa gcggttaatctgctgcgcgataaagggttacttgatgtctggtgacctgggtgattgtcacc cagggcgacgtgatgagtagcgtgggttctactaataaccacgcgtattttaacggtagag taa

Allows you to target upstream regulatory elements

(NT in this case = nucleotide)

Use KEGG database genome map to identify location of your target gene relevant to other genes




Use BLAST (Basic local alignment search tool) to assess off-target binding

Sequence ID: [CP014348.1](#) Length: 4657541 Number of Matches: 163

Range 1: 3891267 to 3891292 [GenBank](#) [Graphics](#)  Next Match  Previous Match


Score	Expect	Identities	Gaps	Strand
52.0 bits(26)	1e-07	26/26(100%)	0/26(0%)	Plus/Plus

Query	1	ATGAAACTCGCCGTTTATAGCACAAA	26
Sbjct	3891267	ATGAAACTCGCCGTTTATAGCACAAA	3891292

Range 2: 392405 to 392417 [GenBank](#) [Graphics](#)  Next Match  Previous Match  First Match

Score	Expect	Identities	Gaps	Strand
26.3 bits(13)	5.8	13/13(100%)	0/13(0%)	Plus/Minus

Query	4	AAACTCGCCGTTT	16
Sbjct	392417	AAACTCGCCGTTT	392405


Range 3: 1595715 to 1595727 [GenBank](#) [Graphics](#)  Next Match  Previous Match  First Match

Score	Expect	Identities	Gaps	Strand
26.3 bits(13)	5.8	13/13(100%)	0/13(0%)	Plus/Minus


Query	1	ATGAAACTCGCCG	13
Sbjct	1595727	ATGAAACTCGCCG	1595715

DNA gel electrophoresis

- Similar concept to SDS-PAGE
- (-) charged DNA migrates through 1% agarose gel toward positive anode
- DNA fragments are separated by size



Diagnostic Digest example



Why are there multiple bands in the uncut plasmid?

Is that OK?

For Today:

- Load and electrophorese DNA gel of diagnostic digest
- Design gRNA primer to improve ethanol or acetate production
 - MUST upload primer to Wiki before leaving class!

For M2D3:

- Read your journal club article and chose 4 figures which are most important to the paper's main conclusion
 - Answer questions on the wiki about those 4 figures
- Email Noreen about presentation day



- Today is a lot of design work on your computer!
- BE Faculty Conference Room is open for you
 - **16-339**
- Office hours begin at 5pm in 16-339
 - Don't need to have specific questions
 - Space to work