

M1D1: Practice cell culture and begin sub-nuclear foci assay

9/11/19

1. Orientation Quiz
2. Prelab discussion
3. Learn about comet chip assay
4. Tissue culture practice
 - a) Seed for gamma-H2AX assay

Office Hours →

Noreen

Wed 10am-12pm
Fri 10am-12pm
in 16-317

Leslie

Wed 9-10am
Fri 4-5pm
in 16-469

Becky

Tues 12-1pm
Thurs 12-1pm
in 16-220



by appointment: nlyell@, lesliemm@, rcmeyer@

Mod 1: Major Assignments

- **Data summary** (15%)
 - In a team, submit on Stellar
 - Draft due 10/14, final revision due 10/26
 - Format: Bullet points, .PPTX
- **Mini-presentation** (5%)
 - Individual, submit video via Gmail
 - Due 10/19 by 10pm
- **Lab quizzes** –be on time!
 - M1D4 and M1D7
- **Notebook** (part of 10% Homework and Notebook)
 - Due 10/4 at 10pm, graded by Colin
- **Blog** (part of 5% Participation)
 - Due 10/15 at 10pm


Elevator pitch

Overview of Module 1: Measuring Genomic Instability

Quantify DNA damage in mammalian cells following exposure to
methylating agent and arsenite


Why these chemicals?


MMS ~~MMN~~




γ H2AX assay: ~~45~~ Stone

- Immunofluorescence
- Cellular response to DNA damage




CometChip assay: 

- single cell gel electrophoresis
- single strand breaks



Tissue culture sterile technique

- **70% ethanol** everything:
 - Wipe cabinet before and after use
 - Wipe everything that enters the cabinet
 - Do not spray cells with EtOH
- **Do not disturb air flow:**
 - Do not block grille or slots
 - Minimize side-to-side arm movements
 - Work > 6" away from sash
 - Leave blower *on always*
- Do not talk into incubator!
- Only open sterile media in hood



Mammalian Cell Culture Medium

We are using CHO cells



Food:

- DMEM (Dulbecco's Modified Eagle's medium)

- **Defined**

sugar, amino acids, vitamins,
salts, phenol red = pH indicator



- FBS (fetal bovine serum)

- **Undefined**

growth factors, lipids
cholesterol



Non-food:

- antibiotics:

- penicillin
- streptomycin

} bacterial growth

Mammalian Cell Culture Terminology

- **Confluence**

Density


- **Splitting**

Subculturing


- **Seeding**

Moving cells to
a new container

Low Density



High Density



80% confluent
↳ split them



General steps for splitting cells

Why do we do each step?

1. Look at cells, estimate confluence
1b. Remove media-aspirator
healthy, confluent media color
2. Rinse with PBS washing off debris / anti-trypsin
3. Detach cells with trypsin (enzyme)
3b. Incubator for 2min
4. Count cells
Know how many
5. "Seed" new culture vessel

Counting cells

10 μ l cell suspension





- Hemocytometer

- Trypan blue \rightarrow cells dead 10%

cells / mL = 10,000 x
average of 4 corners

$$2 \times 4 = 8 / 4 = 2 \times 10,000 = \frac{20,000 \text{ cells}}{\text{mL}}$$



Upcoming assignments: What should go in your notebook?

Laboratory notebook entry component:

Points:

Complete	Partial	Incomplete
----------	---------	------------

Date of experiment (include Module#/Day#) and Title for experiment	1	0.5	0
--	---	-----	---

Hypothesis or goal / purpose	1	0.5	0
------------------------------	---	-----	---

Protocols (link to appropriate wiki sections)	1	0.5	0
---	---	-----	---

Notes on protocol changes / clarifications	1	0.5	0
--	---	-----	---

Observations	2	1	0
--------------	---	---	---

*Visual details			
-----------------	--	--	--

*Qualitative information			
--------------------------	--	--	--

*Raw data			
-----------	--	--	--

Data analysis	3	1.5	0
---------------	---	-----	---

*Calculations			
---------------	--	--	--

*Graphs and Tables			
--------------------	--	--	--

Summary and interpretation of data	3	1.5	0
------------------------------------	---	-----	---

*What did you learn?			
----------------------	--	--	--

*How does this information fit into the larger scope of the project?			
--	--	--	--

Information is clear	2	1	0
----------------------	---	---	---


All days represented	1	0.5	0
----------------------	---	-----	---

OVERALL /15

Due 10pm after each module, as posted on wiki

[http://engineerbiology.org/wiki/20.109\(F18\):_Assignments](http://engineerbiology.org/wiki/20.109(F18):_Assignments)

How should you format your notebook?



The screenshot shows a digital notebook interface with a toolbar at the top containing icons for insert, heading, bold, italic, underline, strikethrough, link, text color, background color, list, indent, and more options. The main content area displays a notebook page with a title, date, and several sections marked with checkmarks.

M1D1: In silico cloning and confirmation digest of protein expression vector

THURSDAY, 2/8

Hypothesis or goal:
What are you testing and what do you expect of your results?

Protocols: [include link to wiki]
Part 2: Construct pRSETb FKBP12 in silico




- Include all work / notes / images / sequences generated.
- Be sure to note any interesting observations or protocol changes!

Part 3: Confirmation digest

- Include completed table with volumes.
- Include calculations.
- Be sure to note any interesting observations or protocol changes!

Summary and interpretations:
What, if any, conclusions can be made and what does this prepare you to do next?

How should you organize your notebook?

- Entitle your project “20.109(F19)_YourName” 
 - Make each module a new folder
 - Make each day a new entry within module folder
- Share the project with Becky and Colin (if you have not already)
 - Right-click and choose ‘settings’ 
 - Add collaborators by email address
 - rcmeyer@mit.edu
 - cyhkim@wi.mit.edu 

Today in Lab

1. TC practice
 - Estimate cell size
 - Split and seed CHO cells
2. Make sure to keep notes in Benchling!
3. Review Cometchip assay and watch Engleward lab JOVE video

<https://www.jove.com/video/50607/cometchip-high-throughput-96-well-platform-for-measuring-dna-damage>



M1D2 HW: Create a template for your benchling notebook and make a M1D2 entry from it. For full credit you must include calculations necessary to complete the table in Part 1 on M1D2.

Condition	EdU (final concentration = 5 μ M)	As (final concentration = 80 μ M)	Media	Final volume (1 mL / well)
No treatment control (no MMS, no As)	1 μ L	0	2 mL	2 mL
MMS, no As				
no MMS, As				
MMS, As				

