

- **Announcements**
- **Pre-lab Lecture**
 - ❖ Your Colony Results
 - ❖ Tissue Culture
 - ❖ Lipofection Workflow
 - ❖ Samples for HR experiment
 - ❖ Today in Lab: M1D6

Announcements

- Slide summary due in one week, Oct 12th
- No class Oct 11th or 12th
- OH change this coming week:
 - Monday 4-5 pm
 - Tuesday... 10-12? 1-3?

Tissue culture tips

- Set up a few inches *behind* the barrier/grate
- Minimize opportunities to bump or expose sterile equipment or your samples
 - Uncap bottles *before* opening pipet
 - Keep tips and dishes *closed* when not in use
 - Avoid passing your hands/arms over open dishes
 - Don't try to hold > 2 things at once! 😊
- Take care not to clog the pipet-aids

Interpreting your ligation results

Group Colour	pCX-EGFP (#)	bkb + lig (#)	bkb + ins, no lig (#)	bkb + ins, lig 1 (#)	bkb + ins, lig 2 (#)
Hypothetical Data	1000	0	2	100	100
Red	1380	5	1	8	2
Orange	1200	10	25	15	150
Yellow	10 000	10	0	80	100
Green	4800	15	3	480	409
Blue	10 ⁶	25	0	90	90
Pink/Purple	8000	1	40	97	188

Consider...

- Why might some groups not have gotten many exptl colonies?
- What does the *no ligase* vs. the *no insert* sample control for? Which one do you expect to have more colonies?
- How big is the variation for identical samples, and what are the possible sources of error causing the variation?

Tissue culture (TC) environment

- What will “feel” physiological to a cell?

pH \sim 7.2-7.4 \leftrightarrow CO₂ (often 5%)

T = 37°C

ambient O₂

humidity

[salts]

sterility*

Tissue culture (TC) medium

- What do cells need to survive?
 - (source for energy - glucose, glutamine, NaPyr)
 - essential amino acids
(non-essential a.a.)
vitamins, minerals, lipids
 - } building blocks
or co-factors
for rxns
 - Serum - cytokines (e.g. growth factor)
 - non-food/optical
antibiotics
pH indicator

Lipofection method

- DNA carrier is similar to the cell membrane
- Efficient transfection (can be >95%)

Figure 6 - Outline of transfection procedure for Lipofectamine™ 2000 Reagent






Figure from Invitrogen website

Lipofection workflow

Wait 5-30 min



... then add to



Wait 20 min



... then add to

Lipofectomine
in Opti-MEM

DNA/ in
Opti-MEM

Lipid/nucleic
acid complexes

$V_{ol} = 50 \mu L$
 $V_L = 2.5 \mu L$


$50 \mu L$ or
 $150 \mu L$

control exp. *triplicate*

16.5x prep,
aliquot

$0.1 \mu g = ? \mu L$


$[DNA] = 0.05 \frac{\mu g}{\mu L}$



Wells with MES cells

Controls for HR assay

- How do you know if your experiment worked?



Experimental samples for HR assay

- How might we increase HR frequency?

• break DNA } before txn - chemical, digest, UV
• vary $\Delta 5:\Delta 3$ } after txn - irradiate, chemical
• vary [DNN] • restriction site location

- Plan for today

– Baseline: $\Delta 5:\Delta 3$ 1:1
– A: $\Delta 3$ cut P_{nucI} - blunt
– B: $\Delta 3$ cut P_{nucII} - sticky
– C: walking on your dishes

Today in Lab: M1D6

- If you were in TC first on D5, go second today
- Half at a time: Lipofection of MES in TC
- Half at a time: Lab practical (**1.5 hours** exactly)
- No need to hand in lab notebooks this week
- Sign up for a Friday FACS time MIDB 'Talk' page
 - you only have to come to lab at your time
 - no lab quiz