#### M1D1:

Learn best practices for mammalian cell culture

- 1. Prelab discussion
- 2. Orientation quiz
- 3. Cell culture exercises




## Mark your calendar!

- Data summary (15%)
  - completed in teams and submitted via Stellar
  - draft due 10/4, final revision due 10/14
- Mini-presentation (5%)
  - completed individually and submitted via Gmail
  - due 10/11
- Laboratory quizzes (collectively 5%)
  - scheduled for M1D4 and M1D7
- Notebook (collectively 5%)
  - one entry will be graded by Aimee 24 hr after M1D7
- Blog (part of 5% Participation)
  - due 10/5 via Blogspot




## Overview of M1: genomic instability

Research question: Does exposure to As inhibit, or decrease, repair of  $H_2O_2$ -induced DNA damage, raising the possibility that combined exposure is an important risk to public health?



#### 1. Use repair foci experiment to measure DNA breaks

 Examine effect of H<sub>2</sub>O<sub>2</sub> +/- As on double strand DNA breaks by measuring γH2AX foci formation



## 2. Use high-throughput genome damage assay to measure DNA damage

 Measure effects of H<sub>2</sub>O<sub>2</sub> +/- As on DNA damage by measuring DNA migration in agarose matrix





## We will use human lymphoblastoid cells

- Specifically, what cell line are we using in M1?
- What are primary cells? Why are they difficult to use in experiments?

Why are cancer cells easier to use in experiments?


 What growth conditions are important when culturing mammalian cells?





## Biosafety cabinets are used to maintain sterility

- Spray everything with 70% ethanol
  - Wipe cabinet before and after use
  - Wipe everything that enters the cabinet
  - Do not spray cells with EtOH
- Do not disturb air flow
  - Do not block grille or slots
  - Minimize side-to-side arm movements
  - Work > 6" away from sash
  - Leave blower on
- Do not talk into cabinet or incubator!
- Only open sterile media in the cabinet



### Growth medium is used to culture cells



#### Food

RPMI 1640 (Roswell Park Memorial Institute)




• FBS (fetal bovine serum)




#### Non-food

Antibiotic solution: penicillin and streptomycin

### Hemocytometers are used to count cells





- Trypan blue mixed with cell suspension at 1:10 ratio, then 10  $\mu$ L added to hemocytometer
- Cells within highlighted sections of the hemocytometer grid are counted


# cells / mL = average # of cells in the 4 highlighted boxes \* 10,000

## The language of cell culture


• Confluence

• Splitting / Sub-culturing

Seeding







High density

### For today...

- Choose a team name!
- Complete Orientation quiz
  - Submit to Stellar by 10 pm
- Work through cell culture exercises
  - Be sure to record your notes in your laboratory notebook

### For M1D2...

- Prepare a template for Benchling laboratory notebook entries
- Be sure to share your Benchling laboratory notebook




## What should go in your notebook?

Laboratory notebook entry component:	Points:		
	Complete	Partial	Missing
Date of experiment (include Module#/Day#) and Title for experiment	1	0.5	0
Mypothesis or goal / purpose	2	1	0
Protocols (link to appropriate wiki sections)	1	0.5	0
Answering questions embedded in wiki sections	5	3	0
Observations from demonstrations and video tutorials	3	2	0
Visual details			
Qualitative information			
Raw data			
Oata analysis	3	2	0
Calculations			
Graphs and Tables			
Summary and interpretation of data	3	2	0
What did you learn?			
How does this information fit into the larger scope of the project?			
Information is clear	2	1	0
All days represented	5	3	0

Be sure to include your responses to the prompts within the laboratory exercises!

## How should you format your notebook?



# M1D1: In silico cloning and confirmation digest of protein expression vector

#### THURSDAY, 2/8

#### Hypothesis or goal:

What are you testing and what do you expect of your results?

#### Protocols: [include link to wiki]

#### Part 2: Construct pRSETb FKBP12 in silico

- Include all work / notes / images / sequences generated.
- Be sure to note any interesting observations or protocol changes!

#### Part 3: Confirmation digest

- Include completed table with volumes.
- Include calculations.
- Be sure to note any interesting observations or protocol changes!

#### Summary and interpretations:

What, if any, conclusions can be made and what does this prepare you to do next?

How should you organize your notebook?

- Title your project "20.109(F20)\_YourName"
  - Make each module a new folder
  - Make each day a new entry within module folder
- Share the project with Instructors and Aimee
  - Right-click and choose 'settings'
  - Add collaborators by email
    - nllyell@mit.edu
    - amoise@mit.edu
    - rcmeyer@mit.edu
    - mebane@mit.edu

