M1D5: Analyze CometChip data and treat cells for sub-nuclear foci assay

09/27/18

- 1. Treat cells for γ H2AX assay
- 2. CometChip data analysis—post on wiki
- 3. Get ahead on reading journal article and homework for next time

Announcements

- Extra office hours coming soon (probably Saturday 10/6)
- Data Summary draft due Monday 10/8 (11 days away!)

Clarifications—Output of Alkaline CometChip Assay



No Damage

- Supercoiled nucleoid
- Little or no migration



High Damage

- SSBs, DSBs, abasic sites, alkali
 labile sites
- forms a "Comet tail"

Genomic damage from direct strand breaks and <u>repair intermediates</u>

Updated slide—Buffers used for CometChips

- Alkaline lysis solution
 - 2.5 M NaCl, 100 mM Na₂EDTA, 10 mM Tris
 - pH 10
 - Triton X-100
- Unwinding & electrophoresis buffer
 - 0.3 M NaOH
 - 1 mM Na₂EDTA
 - pH 13.5
- Neutralization buffer
 - 0.4M Tris
 - pH 7.5
- SYBR Gold DNA stain in PBS



Overview of Mod 1: Measuring Genomic Instability



1. Optimize comet chip assay

Test loading variables



Tod ay

2. Use comet chip assay to measure DNA damage

Measure effects of H₂O₂ on +/- DNA-PK cell lines



- 3. Use immuno-fluorescence assay to measure DNA damage
- Examine effect of H₂O₂ on γH2AX foci formation

Assess DNA damage in tumor cells with & without DNAPKcs

Treat captured cells in comet chip with H_2O_2 (oxidative damage)



Lyse cells & unwind DNA (DNA still captured agarose in overlay)

Agarose Electrophoresis







Stain DNA and image via fluorescence microscopy

Measuring DNA damage via yH2AX Assay



γH2AX = phosphorylated H2AX histone, indicative of DSBs (and potentially other types of DNA damage)

Fix cells and stain with antibody that marks γH2AX



Blue: DNA Green: γH2AX staining

Using immunofluorescence: γH2AX assay to detect double-strand DNA breaks



- Histone H2AX phosphorylated at Ser139 if DSB
- Use antibodies against γH2AX

protein of interest	🔺 γΗ2ΑΧ		
primary antibody	Mouse anti-human anti-γH2AX		
secondary antibody	👗 goat anti-mouse		
Fluorophore (conjugated to secondary antibody) exc./ em. wavelengths	488 / 520 nm		

Treating cells with H_2O_2 for $\gamma H2AX$ assay



Practically using immunofluorescence: γH2AX assay to detect double-strand DNA breaks



What do we hypothesize will happen...



30min or thr

Detecting Repair over time in Human Cells

Each color represents a different human cell line's response to mutagen after an initial exposure followed by recovery time

Oxidative Damage

Alkylation Damage



- Consider difference in kinetics of seeing repair in CometChip vs γH2AX assays
- Will have a paper discussion next week with repair kinetic data

from Prof. Engelward's lecture slides

CometChip Data analysis in ImageJ and MATLAB

1. ImageJ

- from several images per well to one stack per well
- GenImageStacks_single image.txt

2. MATLAB

for each comet in stack, calculates intensity of head and tail, as well as length of tail

3. Excel

- export data from Matlab and compile
- post data to the wiki

T/R purple

buffer_data enzym





14.49

18.78

12.14

10.79

10.36

%Tail DNA

7.33

8.73

11.93

10.54

10.51

%Head DNA

10.03

11.94

10.77

9.76

9.53

8.59

6.86

10.37

14.10

15.28



- Cometnumbers: how many comets were used for calculation in each well (= stack)
- %Head DNA = 100 * HeadFluorescence / (HeadFluorescence + TailFluorescence)
- %Tail DNA = 100 * TailFluorescence / (HeadFluorescence + TailFluorescence)

26.58

34.69

11.85

11.67

OTM (um)

9.68

• Olive tail moment (OTM) = (%TailDNA / 100) * (TailCenterOfMass – HeadCenterOfMass)

triplicates

Tail Len. (um)

Comet Len. (um)

37.04

37.87

10.32

9.29

• Tail length

в

С

D

Ε

F

ometnumbers

Comet length

Major assignments for Mod1

- Data summary draft
 - due by 10pm on Mon., October 8
 - revision due by 10pm on Sat., October 20

Summary content

- 1. Title
- 2. Abstract
- 3. Background & Motivation
- 4. Figures, Results & Interpretation
- 5. Implications & Future Work
- Mini presentation due by 10pm on Sat., October 13
- Blog post for M1 due by 10pm on Tues., October 9

M1 Data summary Architecture



HW M1D6: Implications & Future Works Implications and Future Work: potential topics [edit] Revised methods

- Topic: Did your results match your expectations?
 - · If no, provide a putative explanation. If yes, how can you further test if your hypothesis is correct?
- Topic: Based on the results, whether they matched your expectations or not, what experiments might you recommend next?
 - Follow-up experiments could distinguish between competing explanations of a given outcome or broaden the sample set for a
 question you already asked, to give just two examples.
- Topic: How might this assay be improved?
- Topic: How might this assay be used as a research tool? in the clinic? in industry?

Specific to your research

Check for email from Noreen tomorrow !

In your Data summary tie together (and mirror) your background and motivation, and implications and future work

Broad in scope

HW M1D6: Mini Presentation Outline

- Follow time and content guidelines
- Introduce yourself and your research
- Clearly state your hypothesis to identify main question
- Be quantitative when stating your findings (NOT "This was more/less than...")
- For this HW assignment put placeholder statements for key findings

Category	Approximate worth	Elements of a strong presentation
Content	50%	 Did you introduce your research? Did you include the key findings (and the techniques used to gather these results, if necessary)? Was the importance of your project clear?
Organization	25%	 Is the presentation logical and easy-to-follow? Are the main points emphasized? Did you include transition statements such that the presentation 'flows' and is easily followed/understood?
Delivery	25%	 Do you show confidence and enthusiasm? Did you use appropriate language (technical or informal, as appropriate)? Is your speech clear?

HW M1D6: Prepare for in-class paper discussion

- Consider discussion guidelines on wiki while reading the paper
- Contributing to the discussion is impt. for your participation score

Cell Cycle 12:6, 907-915; March 15, 2013; © 2013 Landes Bioscience

Single-cell microarray enables high-throughput evaluation of DNA double-strand breaks and DNA repair inhibitors

REPORT

David M. Weingeist,^{1,†} Jing Ge,^{1,†} David K. Wood,² James T. Mutamba,¹ Qiuying Huang,³ Elizabeth A. Rowland,¹ Michael B. Yaffe,^{1,3,4,5} Scott Floyd^{4,6} and Bevin P. Engelward^{1,*}

In lab today

- 1. Choose recovery time—check in with me
- 2. Treat cells for γ H2AX and set timer for recovery—fix cells at appropriate time.
- 3. Obtain CometChip data from me and finish analysis today—post results on wiki.
- 4. With extra time, consider reading journal article and doing homework for next time.