

Finish this (start PCR) by 1:25pm

1. Make 100 uM gRNA stock (e.g. 20 nmol → 200 uL nuclease-free water)
2. Make primer mix (10 uM of each primer in 100 uL volume)
 - 10 uL gRNA (forward primer)
 - 10 uL RevP (reverse primer)
 - 80 uL nuclease-free water
3. Assign your gene a color
4. Take PCR tube and label with color you chose
 - 10.25 uL nuclease-free water
 - 1.25 uL primer mix (from #2)
 - 1 uL pgRNA plasmid
 - 12.5 uL Q5 Hot Start 2x Master Mix



M2D3: Generate gRNA plasmid


10/18/19

1. Start PCR
2. BE Communication workshop: Journal Club presentations, 1:30pm, 16-220
3. Pre-lab discussion
4. Set up reaction to generate gRNA_target plasmid
5. Transform pgRNA_target into bacteria


Announcements

- Mini Presentation due Saturday (10/19)! → email to gmail
- Journal club starts on Tuesday!
 - See evaluation rubric on wiki
 - If coming to office hours to prep, email us ahead of time so we can be sure to read your paper

Mod 2 experimental overview



Today's goal: make psgRNA_target



Insert gRNA sequence into expression vector to make pgRNA_target using site directed mutagenesis

Note: sgRNA = gRNA


Use site-directed mutagenesis (SDM) to engineer plasmid DNA

- NEB Q5 SDM kit
- Create specific, targeted changes in double-stranded plasmid DNA

① Forward primer:
*contains mutation
GOI*


• PCR product: *Linear*

• Final product:
circular plasmid



Insertion of DNA via SDM

gRNA sequence (circled) FWD primer pdCas9 tag
Rev. primer







After PCR




After blunt ligation

SDM Part 1: PCR amplification of DNA

Ingredients	
	FWD primer (gRNA target seq)
	REV primer (provided)
Template = <u>pgRNA</u> ?	
dNTPs	
	Polymerase
Buffer (pH, cofactors like Mg^{2+})	
H_2O	




SDM Part 2: Recover circular plasmid product using Kinase-Ligase-Dpn1 (KLD) enzyme mix



 Phosphorylation video:

<https://www.neb.com/tools-and-resources/video-library/the-mechanism-of-dna-phosphorylation>




K: phosphorylate PCR product ends

L: repairs phosphate backbone of plasmid

D: removes methylated DNA (template removal)


We have covered two ways to engineer DNA

“Traditional” plasmid cloning
by restriction enzyme digest



pdCas9

Site directed mutagenesis



pgRNA_target

M2 major assignments

- **Research Article** (20%)
 - Individual, submit on Stellar
 - ~~Word~~ Word document
 - Due 11/12 by 10pm (no revision)
- **Journal Club Presentation** (15%)
 - Individual, during lab section, video recorded, 16-336
 - Watch the video/discuss with Noreen after
 - Powerpoint slides due 1pm on Stellar, on the day of presentation
- Lab quizzes
 - M2D5 and M2D8
- Notebook (part of 10% Homework and Notebook)
- Blog (part of 5% Participation)
 - By 11/1 (Journal club reflections)
 - By 11/12 (Mod2 material/research article)

Tuesday (10/22)

Journal Club I

- Submit presentation slides to Stellar by 1pm
- Presentations should be 10min
 - PLEASE practice your talk out loud at least once!
- Tell us a narrative from the paper
 - Can't present all the data
 - Allowed to pull figures directly from paper to put on slides (title slide = citation)
- You will present from a mac (Keynote or Powerpoint ok)
- Q&A will start with student questions
 - Asking questions counts toward your participation grade
- There will be SNACKS (Trader Joe's– Leslie buys & takes requests)

Friday (10/25)—start homework early

M2D5 HW: Intro, Schematic, Discussion

- Draft Introduction
 - Draft the entire first “Big Picture” paragraph
 - Overview/ topic sentence of each additional paragraph
 - References in text and brief summary of each reference
- Schematic of Mod2 ***theoretical/experimental approach***
 - not experimental overview
 - Create image (do not take and reference published schematics)
 - Include a figure title and caption
- Draft Discussion for confirmation agarose gel figure
 - Draft a paragraph

Reporting and interpreting your data

RESULTS

1. What was the overall goal of these data?
 - State concisely as an introductory sentence.
2. If applicable, what was the result of your control?
 - Was it expected?
3. What was your result?
 - Was it expected?
4. What does this motivate you to do next?
 - Specifically, what experiment follows?

DISCUSSION

1. What evidence do you have that your result is correct or incorrect?
 - How do your controls support your data?
2. In sum, what do your data suggest or indicate?
 - Do your data support your hypothesis? Why?
3. What does this motivate you to do next?
 - Specifically, what is the next research question?

Today in lab...

1. Continue your gRNA insertion/amplification reaction using reagents at front bench
2. Work on M2D5 homework and Journal Club presentations in down time