

M2D2: Design gRNA for CRISPRi


10/18/17

- Pre-lab discussion
- Part 1: Load digests from M2D1 in agarose gels
- Mid-lab discussion: gRNA design considerations
- Part 2: Select target gene in fermentation pathway
- Part 3: Design gRNA for CRISPRi system

From Last Time: Ligation Calculation

Optimal backbone-to-insert ratio

- Ideally, want 4:1 insert:backbone
 - molar ratio, **not** mass or volume




- What if too much insert?

Multiple inserts daisy-chained, inserted into backbone

- What if too much backbone?

Multiple backbones ligating to one another




Today Part 1: DNA: negatively charged

Separate DNA by gel electrophoresis

Agarose gel electrophoresis

- Driving force:
charge, electric field

- Separates DNA by:
size





Today Part 1:

Visualize DNA + save a picture!

- DNA Loading dye (6X):
 - 1) Bromophenol blue
 - small molecule, runs at ~500 bp
 - Visual marker of progress
 - 2) glycerol
 - viscous and dense
 - DNA sinks
- Sybr-Safe DNA stain:
 - DNA intercalator
 - Fluorescent
 - Visualize using UV or blue light

➤ Safety : wear nitrile gloves




Let's do this now...

1. Pick up digest from the front bench and prepare samples to be loaded on the agarose gels
2. Two groups at a time can load their samples and we will start gels after 2 groups per gel are loaded
3. Regroup for another short discussion before moving onto Part 2

Today's goal:


Design gRNA sequence to repress a gene, such that the production of ethanol or acetate will increase.

CRISPRi system overview




- Target gene increase ethanol or acetate production
- pgRNA_target express a gRNA to specific gene
- pdCas9 block expression

CRISPRi system can block the RNA polymerase




How would you use CRISPRi to increase yield of the desired product?




NAD^+

Forster and Gescher. 2014

Let's review binding partners:




Which region of the gene will you target?



Design of gRNA for CRISPRi system

- (1) If you target the template DNA strand, the gRNA sequence will be the same as the transcribed (nontemplate) sequence.
- (2) If you target the nontemplate strand, the gRNA sequence will be the reverse-complement of the transcribed (template) sequence.




Practically: Using KEGG Database

Position	1937649..1939091 Genome map
AA seq	480 aa AA seq DB search MSRRLRRTKIVTTLGPATDRDNNLEKVIAGANVVRMNFSHGSPEDHKMRADKVREIAAK LGRHVAILGDLQGPKIRVSTFKEGVFLNIGDKFLLDANLGKGEVDKEVKVIDYKGLPAD VVPGDILLDDGRVQLKLEVQGMKVFTETVGGPLSNNKGINKLGGGLSAEALTEKDKA DIKTAALIGVVDYLAVSFPRCGEDLNYYARRLARDAGCDAKIVAKVERAEAVCSQDAMDDII LASDVVMVARGLGVEIDGPTELVGIQKALIRRARQLNRavitatQMMESMITNPMPTRA VMDVANAVLDGTDAVMLSAETAAGQYPSETVAAMARVCLGAEKIPSINVSKHLDVQFDN VEEAIAMSAMYAANHLKGVTIAITMTESGRTALMTSRISSGLPIFAMSRHERTLNLITALY RGVTPVHFDSDANGVAAASEAVNLLRDKGYLMMSGDLVIVTQGDVMSTVGSTNTTRILTVE
NT seq	1443 nt NT seq +upstream 0 nt +downstream 0 nt atgtccagaaggcttcgcagaacaaaaatcggttaccacgttagggcccagcaacagatcg gataataatcttggaaaaaatgttatcgccgggtgcacgttgcgtatgaactttct cacggctcgcttgcagaatccacaaaatgcgcgcggataaaagtgcgttagattgcgc ctggggcgcatgtggctattctgggtgcacccgcggggccaaaatccgtgtatccacc tttaaagaaggcaagtttccatattggggataaaattccctgcgtgcacgcacacttgc ggtaaagggtgaaggcgacaagaaaaatgcgttatcgactacaaaggcgcgcgcgc gtcgtgcctgggtgacatcctgcgtgcgtgcgtgcgtgcgtgcgtgc gttcaggggcatgaaagtgttaccgcgactaccgcgcgtggccctctccaacaataaa ggatatcaacaacttggccgcgggttgcgttgcgttgcgttgcgttgcgttgcgttgc gacattaagactgcgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc ggcgaagatctgaactatgcgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc gttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc ctcgccctgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc gaactggtcggcattcagaacggcgatggatccatggatccatggatccatggatcc atcacggcgaccatgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc gtcatggacgttagcaaaccgcgttgcgttgcgttgcgttgcgttgcgttgc actgcgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc gcggaaaaatcccgagcatcaacgcgttgcgttgcgttgcgttgcgttgcgttgc gtggaaaggatattgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc gcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc gttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc gcggtaatctgtcgccgataaaaggtaacttgcgttgcgttgcgttgcgttgc caggcgacgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgcgttgc taa

NT here means nucleotide sequence

Practically: Using KEGG Database

Genome Map



operon: a group of genes transcribed together,
controlled by single promoter

Practically: Using BLAST (Basic local alignment search tool) to assess off-target binding

Sequence ID: [CP014348.1](#) Length: 4657541 Number of Matches: 163

Range 1: 3891267 to 3891292 [GenBank](#) [Graphics](#)

▼ Next Match ▲ Previous Match

Score	Expect	Identities	Gaps	Strand
52.0 bits(26)	1e-07	26/26(100%)	0/26(0%)	Plus/Plus

Query 1 ATGAAACTCGCCGTTTATAGCACAAA 26
Sbjct 3891267 ATGAAACTCGCCGTTTATAGCACAAA 3891292

Range 2: 392405 to 392417 [GenBank](#) [Graphics](#)

▼ Next Match ▲ Previous Match ▲ First Match

Score	Expect	Identities	Gaps	Strand
26.3 bits(13)	5.8	13/13(100%)	0/13(0%)	Plus/Minus

Query 4 AAACTCGCCGTTT 16
Sbjct 392417 AAACTCGCCGTTT 392405

Range 3: 1595715 to 1595727 [GenBank](#) [Graphics](#)

▼ Next Match ▲ Previous Match ▲ First Match

Score	Expect	Identities	Gaps	Strand
26.3 bits(13)	5.8	13/13(100%)	0/13(0%)	Plus/Minus

Query 1 ATGAAACTCGCCG 13
Sbjct 1595727 ATGAAACTCGCCG 1595715

M2D3 HW: Figure/Caption/Results

- Figure = agarose gel image with title and caption
 - don't add drawings/modification on top of image, next to image is sufficient
- Results **paragraph** should follow below figure and results subsection title (take-home message)
- Minimize interpretation in the results section in Mod2 report (separate discussion section; *M2D5HW*)
 - State the result of the experiment without discussing conclusions drawn from experiment
- Figure/text expectations same as Mod1 report, formatting major difference
 - all text in paragraph form

In lecture discussion Thursday

Multiple Gene Repression in Cyanobacteria Using CRISPRi

Lun Yao, Ivana Cengic, Josefina Anfelt, and Elton P. Hudson*

KTH—Royal Institute of Technology, Division of Proteomics and Nanobiotechnology, Science for Life Laboratory, Stockholm SE-171 21 Sweden

S Supporting Information

ABSTRACT: We describe the application of clustered regularly interspaced short palindromic repeats interference (CRISPRi) for gene repression in the model cyanobacterium *Synechocystis* sp. PCC 6803. The nuclease-deficient Cas9 from the type-II CRISPR/Cas of *Streptococcus pyogenes* was used to repress green fluorescent protein (GFP) to negligible levels. CRISPRi was also used to repress formation of carbon storage compounds polyhydroxybutyrate (PHB) and glycogen during nitrogen starvation. As an example of the potential of CRISPRi for basic and applied cyanobacteria research, we simultaneously knocked down 4 putative aldehyde reductases and dehydrogenases at 50–95% repression. This work also demonstrates that tightly repressed promoters allow for inducible and reversible CRISPRi in cyanobacteria.

