

# M2D1: Prepare cells for RNA purification

1. Prelab discussion
2.  $\frac{1}{2}$  class to TC to seed cells for RNA purification
3.  $\frac{1}{2}$  group paper discussion of Dietlein et al.

# Office hours

## Noreen

- M 2-5pm
- in 16-317

## Leslie

- W and F 4-5pm
- in 56-341c or lab

## Josephine

- T 4-5pm
- R 10-11am
- in 56-341c or lab

## + Extra

**Saturday, 03/10**

- **10am-5pm**
- **in 56-302**


**Please email us if you can't make office hours and we will schedule a time to meet!**

# Sign up for journal club


- Pick 1 of 24 papers, or suggest your own
- Present M2D6 (April 4th) or M2D7 (April 6th)
- Sign up by adding your name next to paper [**Josephine Bagnall/WF/Color**]
  - first come first serve!
  - you **cannot** switch paper after M2D4 (March 21<sup>st</sup>)
  - only one T/R presenter and one W/F presenter per article

Slot	Day 6 (T/R)	Day 7 (T/R)	Day 6 (W/F)	Day 7 (W/F)
1	<input type="text"/>			
2				
3				
4				
5				
6				
7				

# M2: Experimental overview



Analyze RNA-seq results




DAY 4: Evaluate altered gene expression


DAY 5: Investigate public databases

DAY 8: Identify regulatory motifs

**Examine effect of drug treatment on cell survival**



# Our cell lines: DLD-1 and BRCA2-/-




- DLD-1 = parental
  - from the colon of a male with colorectal adenocarcinoma
- BRCA2-/- = mutant
  - disruption of exon 11 from BRCA2 gene
  - deficient in DNA repair (by homologous recombination)

*Note:* olaparib is a PARP inhibitor (chemotherapy)

# Tissue culture sterile technique

- **70% ethanol** is your BFF:
  - Wipe cabinet before and after use
  - Wipe everything that enters the cabinet
- Do not disturb air flow:
  - Do not block grille or slots
  - Minimize side-to-side arm movements
  - Work > 6" away from sash
  - Leave blower *on*
- Do not talk into incubator!
- Only open sterile items in hood



# Mammalian cell culture medium

What do cells need to survive? *growth, division, viability*



- *defined* RPMI 1640 (Roswell Park Memorial Institute)
  - (a lot of phosphate)
  - often used to culture lymphoid cells

*· glucose*

*· salts*

*· phenol red (pH indicator)*

*· amino acids*

*· vitamins*

*undefined*

- FBS: fetal bovine serum

*· growth factors*

*· lipids*

*· cytokines*

*· cholesterol*

- *defined* antibiotics:

- penicillin

- streptomycin

*} prevent bacterial growth*



# Mammalian cell culture terminology

- confluence =  
*density*


*at ~80% , split cells*

- splitting  
*Passaging, subculturing*


- seeding

*~20%-40%  
on new dish*

Low Density



High Density



DLD-1 cells



# General steps for splitting cells **+WHY?**

1. Look at cells, estimate confluence

*estimate growth/viability*

2. Rinse with PBS

*wash debris / anti-trypsin agents / serum (FBS)*

3. Detach cells with trypsin


*breaks substrate - cell adhesions*

4. Count cells




*seed specific # in new flask*

5. "Seed" new culture vessel

*room to divide & grow*



# Counting number of cells



- Hemacytometer
- Trypan blue  
*stains dead cells*
- # cells / mL = 10,000 x  
average of 4 corners

$$4 \times 10,000 = 40,000 \text{ cells / mL}$$

# Today in lab:

## 1. Tissue Culture (TC)

- 1<sup>st</sup>: Red, Orange, Green, Blue
- 2<sup>nd</sup>: Pink, Purple, Silver
- Protocols printed for TC use, no need to move laptops etc.
- Do not wear PPE in or out of TC room

## 2. Paper discussion of Dietlein *et al.*

- Homework due M2D2
  - Sign up for journal club day (sign up for article next week)
  - Create a single slide from Dietlein *et al.*
- Don't forget about Mod1 assignments!

# M2D2 HW: JC Slide

- Slide= Standard 4:3 powerpoint slide
- Title has a message (not just the figure title)
- Don't put too much on one slide, (1 slide=1 message)
- Don't fill slide with text
- Don't include the caption from paper or a citation