

M2D1: Complete *in silico* cloning of dCas9

10/12/16

1. Pre-lab Discussion
2. Design primers to dCas9
3. *In silico* PCR amplification, digest, and ligation
4. REAL diagnostic digest of pdCas9

(Almost) done with Mod1!

- Mini-presentation
 - due 10pm on Saturday, October 14
 - You don't have to be exhaustive;
tell a focused story
- Data summary
 - ✓ draft due 10pm on Wednesday, October 11
 - receive all comments by Monday, October 16th
 - revision due 10pm on Sunday, October 22nd
- Blog post
 - due 10pm on Monday, October 23



Sign up for journal club (M2D2 homework)

- Pick 1 of 20 papers, or suggest your own
- Present M2D4 (October 24) or M2D6 (October 31)
- Sign up by adding your name next to paper [LM/TR/Rainbow]
 - first come first serve!
 - only one T/R and one W/F student per article
 - Don't pick a paper randomly


Slot	Day 4 (T/R)	Day 6 (T/R)	Day 4 (W/F)	Day 6 (W/F)
1				
2				
3				
4				
5				
6				

“insert”
“fragment”
dCas9


“vector”
“backbone”
p_LtetO-1

How is DNA engineered?

- PCR amplification of DNA:
 - polymerase chain reaction
 - specific primers to the gene of interest
- digestion:
 - restriction enzymes
 - endonucleases
 - sticky ends (or blunt ends)
- ligation:
 - ligase: seal backbone, phosphodiester bond




How do you design primers? **all primers 5' to 3'**




Primer design guidelines

- Length
 - 17-28 base pairs
 - **long enough to be specific**
 - **short enough so T_m not so high**
- GC content
 - 40-60%
 - GC clamp at ends
 - **AT=2, GC=3; GC more stable**
- T_m(primer)
 - < 65 °C
- Avoid Secondary structures
 - hairpins
 - complementation w/i primer sequence
- Avoid repetitive sequences
 - di-nucleotides > 4 (ex. ATATAT)
 - runs > 4 bp (ex. GATGGGGG)



Using PCR to generate *dCas9* flanked by restriction enzyme recognition sites

3 major steps in the PCR cycle: which temperature and why?

- Melt
 - 95°C
 - dsDNA denatures
 - Anneal
 - $T_m(\text{primer}) = 1/2 \text{ primers annealed to target}$
 - $T_{\text{anneal}} \sim T_m(\text{primer}) - 5^\circ\text{C}$
 - Extend
 - 72°C (for Taq)
 - 1 min / 1000 bases
- 

The diagram illustrates the PCR process. At the top, a circular blue line represents the dCas9 gene, with two restriction enzyme sites labeled 'BglII' and 'XhoI' indicated by arrows. A large grey arrow points down to a linear blue line representing the amplified DNA product, which also shows the 'BglII' and 'XhoI' sites flanking the dCas9 gene.


 - Template DNA
 - dCas9 gene
 - Primers
 - specific match to dCas9
 - option to add bases, *e.g.* endonuclease recognition sequence

My favorite PCR animation..


learn.genetics.utah.edu



Digestion: dCas9 insert and plasmid backbone





Ligation: dCas9 insert and plasmid backbone after digestion



Confirmation digest considerations

- Do you have access to the enzymes?
- Are the two enzymes compatible?
- Are fragments easily distinguished on an agarose gel?



Today in lab

1. Reproduce *in silico* (in Benchling) the cloning of pdCas9
 - design primers that would amplify the gene dCas9
 - depict PCR amplification product
 - digestion of dCAS9 PCR product and vector by restriction enzymes
 - ligation of insert and vector ** at 2:30pm we will all work through this calculation together
2. Set up confirmation digests of pdCas9 for agarose gel electrophoresis
 - choose restriction enzymes for diagnostic digest
 - calculate volumes of digest components
 - “set-up” digest and leave overnight at 37°C