M1D7: Complete sub-nuclear foci assay staining and data analysis


10/4/18

- 1. Stats
- 2. Finish staining γH2AX foci
- 3. Visit microscope
- 4. γH2AX image analysis
- 5. Data analysis—ask questions while you're here!

<u>Announcements</u>

- Notebook due 10/5, 10pm
 - Graded in detail: M1D5
- Regular office hours:
 - Fri. (10/5) 2-3pm, 56-341c
- Extra office hours:
 - Sat. (10/6) 10am-12pm, 56-302
 - Mon. (10/8) 12-5pm, 56-302
 - Prof Engelward: Wed. (10/10),
 9am-12pm, 1-3pm
- Data Summary draft due 10pm, Wed. 10/10

Practically using immunofluorescence: γH2AX assay to detect double-strand DNA breaks




Complete vH2AX assay today


Phalloidin stain (F actin)


Together with Secondary Ab

DAPI stain (DNA)

In mounting medium








Mount on glass slide


γH2AX Data Analysis

What is the total amount of vH2AX signal of each nucleus?

- Use DAPI channel to identify nude
- Use ImageJ to quantify total FITC (<u>xtex</u>) fluorescence in each nucleus
- Imaged using same exposure times per channel across all










γH2AX Data Analysis using ImageJ

- Pick a threhold value in the DAPI channel to identify nuclei—typically good to be consistent, use same threshold on all images.
- Identify nucleus location using PAPL channel, & quantify fluorescence in <u>fitc</u> channel.
- Normalize fluorescence intensity by <u>avea</u> of nucleus





X post on wiki Perform analysis on resultant excel spreadsheet



Solidity

0.951

0.979

0.97 0.978

0.96

0.837

Round

0.723

0.929

0.658

0.721

0.485

0.411

1.382

1.077

1.519

1.386

497394

6064824

3901340

5790485

#6	much region	Fred	E	(J.	, vame		2	
	A		В			4	D	
1		Label	1		Area	$^{\prime}$	Mean	
2	1	J_H2O2	_001-00	002	$\overline{}$	10213	487.02	21

2 J_H2O2_001-0002

3 J_H2O2_001-0002

4 J H2O2 001-0002

W Z		factoring in							
, et			Kage Witon), ,					
D	E	F	G	Н					
Mean	Min	Max	Circ.	IntDen	RawIntDen				

1062

1005

549

865

5 J_H2O2_001-0002 10959 659.783 439 2451 0.729 7230565 7230565 2.062 6 J_H2O2_001-0002 38645 595.773 398 1937 0.456 23023658 23023658 2.432	Compare:					· Cometchip VsyHZAX						
F 1 11000 001 0000 10000 000 000 000 000												

371

391

381

380

Killetics

0.768

0.862

0.811

0.838

4973943

6064824

3901340

5790485

i.

- Cell lives
- A202 treatment

10249

8482

11661

591.748

459.955


496.568


arbitrary units


avea


Make stacked images

- Use ImageJ to overlay and color the images
- Don't forget a scale bar: 0.17 um/px









Actin (TxRed)

Nuclei (DAPI)

YH2AX (FITC)

Major assignments for Mod1

- **Data summary** draft
 - due by 10pm on Wed., October 10
 - revision due by 10pm on Sat., October 20

Summary content


- 1. Title
- 2. Abstract
- 3. Background & Motivation
- 4. Figures, Results & Interpretation
- 5. Implications & Future Work
- Mini presentation due by 10pm on Sat., October 13
- Lab notebook for Mod 1 due by 10pm on Fri. October 5
- Blog post for M1 due by 10pm on Thurs., October 11

20.109 Blog Post


- You will receive an invitation to join the class blog
- Possible topics listed on the blog
- Details about use:
 - Do not publish MIT logo
 - Do not post photographs with names tagged
 - Do not write malicious comments
 - Do not plagiarize







M1 Data Summary



Title: take-home message Abstract: the only page *not* in bullet points

ALL bullet points:

- -background and motivation (include references)
- schematics
- -Results and interpretation
- Cell loading
- Comet Chip analysis (team expt.)
- Comet Chip repair data (Engelward lab)
- H2AX assay (team & class data, images)
- Schematics?

Implications and future work


References (see wiki for format suggestions)

Background & Motivation

- Impact statement
 - general background
 - describe previous work in the field
- Specific background (e.g. BER, DNAPKcs, CometChip)
 - introduce topics, pathways and specific technologies necessary to understand the experiment
 - narrow focus to the specific question addressed in your study

Knowledge gap/statement of problem

- what is unknown, therefore motivating your study
- · Hypothesise Be specific to your experiments
 - what do you propose will be the outcome of your study
- A brief preview of your findings
 - Here we show...
 - end with broad implications of the study



The meat of your paper: Results & Interpretation

- Figures and captions
 - Decide on these first
 - Use subpanels
 - Text: limited on figure, explicit in caption
 - reasonable size
 - descriptive title
 - intro sentence in caption
 - caption descriptive of image, very light on methods
- Result bullets
 - Goal / intent / purpose of experiment = intro topic sentence
 - What you did: experiments and expectations, including controls
 - What you found: quantitatively describe your result, referring to the figure
 - What does this indicate: interpret your result, what does it mean?
 - What does this motivate you to do next: transition to next experiment

Example Results slide (from Wiki)




Figure 1: Development of BRET assay to monitor EGFR and SH2 domain interactions. CHO-K1 cells were transfected with Citine-EGFR (A) and renills lucferase (RLuc)-tagged SH2 domains from P₂Cp, Gr62, CTEN, and SHc3 (B). Western bits of CHO-K1 hysters were probed with anti-EGFR (A) or anti-RLuc (B) antibodies. Arrowheads indicate the expected molecular weight of the RLuc-SH2-GR2, CR2 proteins; (1) RLuc-SH2-GF2, CR2 (CR2 RLuc-SH2-GF2, CR2 RLuc

BRET system effectively measures EGFR activation:

- To determine if the BRET system could be used to monitor EGFR activation, CHO-K1 cells were transfected with fluorescent EGFR and luciferase-tagged SH2 domains and a BRET assay was performed after growth factor stimulation.
- CHO-K1 were transfected with Citrine-EGFR in all conditions as indicated by correct molecular weight band at 150 kDs (Figure 1A).
- Several protein bands are present in Mock transfection lane suggesting off-target binding of the RLuc antibody (Figure 1B).
- RLuc alone, RLuc-SH2-Grb2, and RLuc-SH2-CTEN were successfully transfected as indicated by correct molecular weight bands (Figure 1B).
- RLuc-SH2-PLCg and RLuc-SH2-Shc3 did not appear by Western blot analysis bands different from those in the Mock lane are not identifiable. This outcome could be due to protein expression levels below the detection limit by Western blot or to unsuccessful transfection of cDNA.
- BRET signal increased in cells transfected with Citrine-EGFR and RLuc-SH2-CTEN versus Citrine-EGFR and RLuc alone after EGF stimulation. This difference suggests that the BRET signal is specific for an SH2-EGFR interaction versus randomly localized RLuc.
- In sum, these data suggest that the RLuc-SH2 constructs can be utilized to monitor EGFR phosphorylation, as SH2 domain-EGFR association occurs only at sites of EGFR tyrosine phosphorylation. Next, we determined the dynamic range of the BRET assay.

- Start with a very similar paragraph to the last paragraph in your Background/Motivation (restate major results and broad implications)
- Follow same order as in Figures/Results
 - Describe your conclusions from your data
 - Describe caveats and suggest remedy
 - Identify unknowns and speculate within reason
 - Don't make huge generalizations or overreach
- Propose future experiments, identify new questions that arise
- Come back to (the same) big picture topic introduced in background

In lab today

- 1. Finish staining γH2AX coverslips
- 2. Stats exercise
- 3. γH2AX image analysis
- 4. Finish data analysis
- 5. Continue working on Data Summary