

M1D5: Analyze CometChip data and treat cells for sub-nuclear foci assay

09/28/18


1. Treat cells for γ H2AX assay
2. CometChip data analysis—post on wiki
3. Get ahead on reading journal article and homework for next time



Announcements


- Extra office hours coming soon (probably Saturday 10/6)
- Data Summary draft due Monday 10/8 (11 days away!)

Clarifications—Output of Alkaline CometChip Assay



No Damage

- Supercoiled nucleoid
- Little or no migration




High Damage

- SSBs, DSBs, abasic sites, alkali labile sites
- forms a "Comet tail"


**Genomic damage from direct strand breaks
and repair intermediates**

Overview of Mod 1: Measuring Genomic Instability




1. Optimize comet chip assay

- Test loading variables



2. Use comet chip assay to measure DNA damage

- Measure effects of H_2O_2 on +/- DNA-PK cell lines




3. Use immuno-fluorescence assay to measure DNA damage

- Examine effect of H_2O_2 on γ H2AX foci formation





Assess DNA damage in tumor cells with & without DNAPKcs

Treat captured cells in comet chip
with H_2O_2 (oxidative damage)




Lyse cells & unwind DNA
(DNA still captured
agarose in overlay)

Agarose Electrophoresis



Stain DNA and image via
fluorescence microscopy



Measuring DNA damage via γ H2AX Assay



no treatment, H_2O_2 treatment, recovery condition

γ H2AX = phosphorylated H2AX histone, indicative of DSBs (and potentially other types of DNA damage)

Fix cells and stain with antibody that marks γ H2AX




Blue: DNA





Green: γ H2AX staining

Using immunofluorescence: γ H2AX assay to detect double-strand DNA breaks




- Histone H2AX phosphorylated at Ser139 if DSB
- Use antibodies against γ H2AX




protein of interest (<i>antigen</i>)	 γ H2AX
primary antibody	 mouse anti-human anti- γ H2AX
secondary antibody	 goat anti-mouse
Fluorophore (conjugated to secondary antibody) exc./ em. wavelengths	 488 / 520 nm

Treating cells with H_2O_2 for $\gamma H2AX$ assay



Practically using immunofluorescence: γ H2AX assay to detect double-strand DNA breaks



What do we hypothesize will happen...

γ H2AX signal

	+DNAPKcs	-DNAPKcs
A	NONE	NONE
B	LOW	Slightly higher than K
C	higher than row B	Same as K / slightly higher

No H_2O_2

H_2O_2 for 20min

Fixed immediately

H_2O_2 for 20min

Recovery in fresh media

Consider:
 * repair kinetics
 * kinetics of γ H2AX
 phos + dephos
 * slightly diff cell conditions
 * 1mM H_2O_2

30min
1hr


② higher than C
 ① slightly lower than C @ 30min
 1hr lower

J slightly more signal than K


Detecting Repair over time in Human Cells

Each color represents a different human cell line's response to mutagen after an initial exposure followed by recovery time

Oxidative Damage



Alkylation Damage



- Consider difference in kinetics of seeing repair in CometChip vs γ H2AX assays
- Will have a paper discussion next week with repair kinetic data

from Prof. Engelward's lecture slides

CometChip Data analysis in ImageJ and MATLAB

1. ImageJ


- from several images per well to one stack per well
- GenImageStacks_single_image.txt

TK6_0H2O2_01A_0020.tif
TK6_0H2O2_01A_0021.tif
TK6_0H2O2_01A_0022.tif
TK6_0H2O2_01A_0023.tif
TK6_0H2O2_01A_0024.tif
TK6_0H2O2_01B_0020.tif
TK6_0H2O2_01C_0019.tif
TK6_0H2O2_01D_0022.tif
TK6_0H2O2_01E_0023.tif
TK6_0H2O2_01F_0024.tif

TK6_0H2O2_01A_.tif
TK6_0H2O2_01B_.tif
TK6_0H2O2_01C_.tif
TK6_0H2O2_01D_.tif
TK6_0H2O2_01E_.tif
TK6_0H2O2_01F_.tif

2. MATLAB

- for each comet in stack, calculates intensity of head and tail, as well as length of tail



3. Excel

- export data from Matlab and compile
- post data to the wiki


160913TK6_H2O2_6Do	
Team	Analyzed Data
T/R red	buffer_data enzym
T/R orange	buffer_data enzym
T/R yellow	buffer_data enzym
T/R green	buffer_data enzym
T/R blue	buffer_data enzym
T/R pink	buffer_data enzym
T/R purple	buffer_data enzym

What's in the final Excel file?

	01	02	03	04	05	06
A	7.45	7.68	11.33	16.49	34.06	29.43
B	8.59	7.33	10.03	14.49	26.58	37.04
C	6.86	8.73	11.94	18.78	34.69	37.87
D	10.37	11.93	10.77	12.14	9.68	11.71
E	14.10	10.54	9.76	10.79	11.85	10.32
F	15.28	10.51	9.53	10.36	11.67	9.29

triplicates

Cometnumbers	%Head DNA	%Tail DNA	OTM (um)	Tail Len. (um)	Comet Len. (um)	+
--------------	-----------	-----------	----------	----------------	-----------------	---



- Cometnumbers: how many comets were used for calculation in each well (= stack)
- %Head DNA = $100 * \text{HeadFluorescence} / (\text{HeadFluorescence} + \text{TailFluorescence})$
- %Tail DNA = $100 * \text{TailFluorescence} / (\text{HeadFluorescence} + \text{TailFluorescence})$
- Olive tail moment (OTM) = $(\% \text{TailDNA} / 100) * (\text{TailCenterOfMass} - \text{HeadCenterOfMass})$
- Tail length
- Comet length

Major assignments for Mod1


- Data summary draft
 - due by 10pm on Mon., October 8
 - revision due by 10pm on Sat., October 20

Summary content

1. Title
2. Abstract
3. Background & Motivation
4. Figures, Results & Interpretation
5. Implications & Future Work

- Mini presentation due by 10pm on Sat., October 13
- Blog post for M1 due by 10pm on Tues., October 9

M1 Data summary Architecture



HW M1D6: Revise methods (*with partner*)

Mini Presentation **Outline** (*individual*)

- **Follow time and content guidelines**
- Introduce yourself and your research
- Clearly state your hypothesis to identify main question(s)
- Be quantitative when stating your findings (NOT “This was more/less than...”)
- For this HW assignment put placeholder statements for key findings

Category	Approximate worth	Elements of a strong presentation
Content	50%	<ul style="list-style-type: none">• Did you introduce your research?• Did you include the key findings (and the techniques used to gather these results, if necessary)?• Was the importance of your project clear?
Organization	25%	<ul style="list-style-type: none">• Is the presentation logical and easy-to-follow?• Are the main points emphasized?• Did you include transition statements such that the presentation 'flows' and is easily followed/understood?
Delivery	25%	<ul style="list-style-type: none">• Do you show confidence and enthusiasm?• Did you use appropriate language (technical or informal, as appropriate)?• Is your speech clear?

HW M1D6: Prepare for in-class paper discussion

- Consider discussion guidelines on wiki while reading the paper
- Contributing to the discussion is impt. for your participation score

Cell Cycle 12:6, 907–915; March 15, 2013; © 2013 Landes Bioscience

REPORT

Single-cell microarray enables high-throughput evaluation of DNA double-strand breaks and DNA repair inhibitors

David M. Weingeist,^{1,†} Jing Ge,^{1,†} David K. Wood,² James T. Mutamba,¹ Qiuying Huang,³ Elizabeth A. Rowland,¹ Michael B. Yaffe,^{1,3,4,5} Scott Floyd^{4,6} and Bevin P. Engelward^{1,*}

In lab today

yellow 30m
green 30m
blue 1hr
~~pink~~ 30m

1. Choose recovery time—check in with me
2. Treat cells for γ H2AX and set timer for recovery—fix cells at appropriate time.
3. Obtain CometChip data from me and finish analysis today—**post results on wiki.**
4. With extra time, consider reading journal article and doing homework for next time.