

M2D7: Induce CRISPRi system

1. BE Communication workshop
2. Pre-lab
3. Analyze sequencing results
4. Prep for dCas9 induction and mixed-acid fermentation

Major assignments for M2

- **Research Article**
 - Due by 10pm on Mon., November 12th

Research Article content

1. Title
2. Abstract
3. Introduction
4. Materials and Methods
5. Figures and Results
6. Discussion
7. References

- Blog post for journal club due 10pm, 11/2
- Blog post for Mod 2 due 10pm, 11/13

Extra office hours


11/10 (Sat):

11/11 (Sun):

Regular office hours


- Noreen: Mon. 2-5pm (16-317)
- Leslie: Th 2-3pm, Fr 12-1pm (56-341c)
- Josephine: W 12-1pm, Fr 2-3pm (56-341c)
- Email us to schedule a different time

M2 experimental overview




Note: sgRNA = gRNA

Using **CRISPRi** to manipulate the *E. coli* fermentation pathway




Design of gRNA for CRISPRi system



Design of gRNA for CRISPRi system

- (1) If you target the template DNA strand, the gRNA (DNA) sequence will be the same as the transcribed (nontemplate) sequence.
- (2) If you target the nontemplate strand, the gRNA (DNA) sequence will be the reverse-complement of the transcribed (template) sequence.



Please add your targeting info to the wiki today


On the [Class Data](#) page

| Team | Ethanol (E) or Acetate (A) | Gene targeted by CRISPRi gRNA | gRNA sequence end) |
|--------|----------------------------|-------------------------------|--------------------|
| red | Ethanol | ack (indirectly, pta) | GTTTTTTTAGCC |
| orange | Ethanol | ldhA | ATTCAACATCAC |
| blue | Ethanol | ackA | TTTTTAGCCACG |

...

| | Locus targeted (eg. beginning of gene, putative promoter, -35 region) | Target template or nontemplate strand |
|----|---|---------------------------------------|
| | | |
| TT | | |
| | | |

M2D3: Generated pgRNA_target by SDM



← CRISPRi universal *amplification* reverse primer

forward primer including crRNA to be inserted ()
dCas9 handle (→)

Analyzing Sequence Information

- Reverse and complement your reverse primer sequence before alignment
- Check whether your target sequence was successfully incorporated into the pgRNA_target plasmid


```

1>---gaattctaagaatctttgacagctagctcagtcctaggtataatactagct-----gttttagagctagaaatagcaag>73
73>---GAATTCATAAGATCTTTGACAGCTAGCTCAGTCCTAGGTATATACCTAGTAAATCCACTTAAGAGGTAAGGTGTGTTTTAGAGCTAGAAATAGCAAG>269
01>CTGGAATTCATAAGATCTTTGACAGCTAGCTCAGTCCTAGGTATATACCTAGTAAATCCACTTAAGAGGTAAGGTGTGTTTTAGAGCTAGAAATAGCAAG>500
1>---gaattctaagaatctttgacagctagctcagtcctaggtataatactagct-----gttttagagctagaaatagcaag>73


74>ttaaaaataaggctagtcggttatcaacttgaaaaagtggcaccgagtcgggtgcttttttgaagcttgggcccgaacaaaaaactcatctcagaagaggat>173
70>TTAAATTAAGGCTAGTCGGTTATCACTTGAAAGTGGCACCAGTCGGTGCTTTTTTTGAGCTTGGGCCCGAACAAAACCTCATCTCAGAGAGAGGAT>369
01>TTAAATTAAGGCTAGTCGGTTATCACTTGAAAGTGGCACCAGTCGGTGCTTTTTTTGAGCTTGGGCCCGAACAAAACCTCATCTCAGAGAGAGGAT>600
74>ttaaaaataaggctagtcggttatcaacttgaaaaagtggcaccgagtcgggtgcttttttgaagcttgggcccgaacaaaaaactcatctcagaagaggat>173

```


- Sanger sequencing traces are also on wiki for your reference



Induction of CRISPRi system with aTc



Anhydrotetracycline (aTc)



E. coli MG1655 genome


- Expressed constitutively:

gRNA

- Expression induced with aTc:

ΔCas9


CRISPRi 'inactive' in absence of inducer




Without aTc

- Only gRNA present
- No (or little) dCas9


aTc induction of pdCas9



- Tet promoter regulates expression of dCas9 gene



CRISPRi 'blocks' gene expression in presence of inducer



dCas9 protein
associates with
gRNA/target gene
to repress target
gene expression

Set up liquid cultures for mixed-acid fermentation and pdCas9 induction

- What are the necessary components?

- bacteria (MG1655 +/- CRISPRi)

- LB

- atc (induce dCas9 expression)

- Screw Cap tubes (anaerobic)

- antibiotics (amp, cam)

- What control conditions will we have?



- MG1655 +/- CRISPRi

- aerobic vs anaerobic

- atc vs. no atc

Set up liquid cultures for mixed-acid fermentation and pdCas9 induction

- Where do we expect most ethanol if hypothesis confirmed?



M2D8 Assignments

- Quiz on M2D8
- Peer-review methods
 - Do not leave today before receiving Methods to peer-review (soft or hard copy?)
 - Either print it out and manually indicate which part each comment refers to, type out comments, and scan to submit online, or hand in hard copy.
 - Or make digital comments with related highlights in Adobe Acrobat

Today in lab...

- Upload gRNA design information to wiki
- Download your sequencing data from discussion tab and align (using ApE software)
- Prepare media for mixed-acid fermentation inoculations