

Module 1 overview

lecture

1. Introduction to the module
2. Rational protein design
3. Fluorescence and sensors

lab

1. Start-up protein eng.
2. Site-directed mutagenesis
3. DNA amplification

PRESIDENT'S DAY

4. Review & gene analysis
5. Protein expression
6. Purification and protein analysis
7. Binding & affinity measurements
8. High throughput engineering



4. Prepare expression system
5. Gene analysis & induction
6. Characterize expression
7. Assay protein behavior
8. Data analysis

Lecture 7: Binding and affinity measurements


- I. Titration analysis
 - A. Estimating K_d & EC_{50} from fluorescence data
 - B. Multisite binding and cooperativity

- II. Techniques for studying binding

Q: Our design task is to alter calcium sensitivity of pericam derivatives. We chose mutants to make, but **how will we determine the results** of our perturbations?



What is going on during a titration measurement?



$$[P]_{tot} = [P_u] + [P_b]$$

$$K_d = ([P_u][Ca^{2+}]^4)/[P_b]$$

We want to measure K_d from fluorescence. P_u and P_b are fluorescent to different degrees, so that the total fluorescence is equal to the sum of contributions from the two species:

$$F = F_u [P_u] + F_b [P_b]$$

Fluorescence measurements can be converted to fractional saturation:

$$\frac{F_{max} - F}{F_{max} - F_{min}} = \frac{\Delta F}{\Delta F_{max}} = \frac{[P_b]}{[P]_{tot}}$$

$\Delta F/\Delta F_{max}$ is equal to the fractional saturation of pericam calcium binding sites, often abbreviated Θ . Θ can be expressed as a function of calcium concentration as follows:

$$\begin{aligned}\Theta &= \frac{[P_b]}{[P]_{tot}} = \frac{[P_b]}{[P_b]K_d/[Ca^{2+}]^4 + [P_b]} \\ &= \frac{[Ca^{2+}]^4}{K_d + [Ca^{2+}]^4}\end{aligned}$$


$EC_{50} = (K_d)^{1/4}$ is the calcium concentration at which half-maximal binding and fluorescence change occurs.

Calcium sensing could realistically be performed for calcium concentrations near the EC_{50} , and our design goal is equivalent to shifting the EC_{50} of the pericam derivatives.

If we had fluorescence data over a range of calcium concentrations, and assuming these equations accurately describe calcium binding to pericam, we could determine the EC_{50} and K_d using a variety of methods:

2. Hill analysis

- substitute $\Theta = \Delta F / \Delta F_{max}$
- set $y = \log[\Theta / (1 - \Theta)]$
- plot y vs. $x = \log([Ca^{2+}])$ for transition region
- slope tells number of cooperative binding sites, “Hill coefficient”
- x intercept is $\log(EC_{50})$



Hill coefficient reflects **cooperativity**, a phenomenon of binding to multiple binding sites on an individual target: **positive cooperativity** means that binding to one site promotes binding to other sites; **negative cooperativity** means that binding to one site depresses binding to other sites.

In this case, the equation that describes a titration curve (fraction of bound sites vs. calcium concentration) is:


$$\Theta = \frac{[B_b]}{[B]_{tot}} = \frac{[B_b]}{[B_b]K_d/[Ca^{2+}] + [B_b]}$$

$$= \frac{[Ca^{2+}]}{K_d + [Ca^{2+}]}$$

Hill coeff. = 1


Simulated binding curves for **BAPTA** and **pericam** compare as follows:

saturation vs. $\log([Ca^{2+}])$




pericam: sharper transition
BAPTA: broader transition

saturation vs. $[Ca^{2+}]$




pericam: sigmoidal transition
BAPTA: smooth transition

Hill plot




pericam: greater slope (Hill coefficient)
BAPTA: lesser slope

The titration behavior we've been discussing for pericam is idealized, because we have been assuming **complete cooperativity** among the calcium binding sites (all four Ca^{2+} ions bind at once). In fact, binding to individual sites *can occur independently*, but each site's apparent affinity depends on whether the other sites are occupied.



No titration curve provides enough detail (features) to accurately fit all the independent equilibrium constants in the “full picture.” One possibility is to settle for a compromise that fits data from CaM-based sensors:




$$\frac{\Delta F}{\Delta F_{max}} = f_1 \times \frac{[Ca^{2+}]^{n_1}}{K_{d1} + [Ca^{2+}]^{n_1}} + f_2 \times \frac{[Ca^{2+}]^{n_2}}{K_{d2} + [Ca^{2+}]^{n_2}}$$

n_1 and n_2 are the Hill coefficients associated with the transitions from P_u to P_i and from P_i to P_b, respectively. f_1 and f_2 are the fractional fluorescence changes associated with these two steps.

Why is modeling pericam with **two calcium-dependent** transitions particularly appropriate?

One- and two-step transition behavior of CaMeleons
(Miyawaki *et al.*):



- your mutations in pericam may affect the apparent K_d s (or EC_{50} s) and Hill coefficients for two transition steps
- the “resolution” of the two steps may be affected—note that the parent pericam appears to have a single transition, in part due to the E104Q mutation (*cf.* curve with black circles above)
- in some cases, your mutations may affect the relative fluorescence changes of the two transition steps

Note on calcium concentrations

13

How is a range of known calcium concentrations be produced?

- note that EC₅₀s are typically around 1 μM; the most useful calcium concentrations will be **near the EC₅₀**
- you may need > 1 μM protein to make robust measurements
- need to make sure that we know the concentration of **unbound calcium**, as opposed to total calcium concentration

$$\frac{\Delta F}{\Delta F_{max}} = f_1 \times \frac{[Ca^{2+}]^{n_1}}{K_{d1} + [Ca^{2+}]^{n_1}} + f_2 \times \frac{[Ca^{2+}]^{n_2}}{K_{d2} + [Ca^{2+}]^{n_2}}$$

We will solve this problem by using a **calcium buffer**

- analogous to a pH buffer, with pCa²⁺ determined by affinity of the buffer (in our case EGTA) and absolute amount of Ca²⁺ present


$$[Ca^{2+}] = \frac{[Ca^{2+} - EGTA]}{[EGTA] \cdot K_a}$$

- with 10 mM Ca²⁺-EGTA/EGTA mixtures, the amount of pericam present is unlikely to have much effect

Binding may be quantified using methods other than fluorescence

Other techniques for titration curve measurement:


- absorbance spectroscopy,
e.g. O₂ binding to Hb



www.oximetry.org/pulseox/principles.htm


This and following
slides not covered,
Just reference

- circular dichroism, *e.g.* Ca²⁺
binding to TrC fragments




Greenfield (2007) *Nat. Protoc.* 1: 2733-41

- surface plasmon resonance, *e.g.* antibody binding to a target




De Santis *et al.* (2003)
Br. J. Cancer 88: 996-1003

- calorimetry, *e.g.* peptide-protein interaction



Schütz *et al.* (2006) *EMBO J.* 25: 4245-52

- analytical ultracentrifugation, *e.g.* ligand-protein binding



Arkin & Lear (2001) *Anal. Biochem.* 29: 98-107