- Announcements
- Pre-lab Lecture
 - Circuit analogy (enrichment)
 - Where we are/going
 - Today in Lab (M2D5)

Announcements

- FNW (and today's assay!) comments

 - when is time is an issue? short +/- Isec = larger to and
 A420 vs. OD600 out of range only measurements matter

 Journal club coming up


 In one week

 in one week

 when is time is an issue? short +/- Isec = larger to and
 isec = larger to and
- Journal club coming up

 - extra OH besides Tuesday? ⇒ \angle \alpha\ldots\ldots
 - poll by Doodle for one-on-one meetings?
- Start thinking about best day for M2 extra OH

Revisiting circuit analogue: main equation



Revisiting circuit analogue: changing R

 $V_out/I_in = -R$

Shine light on photodiode →

I_in flows →

V_out decreases (- sign) →

Bigger voltage drop →

LED lights up!


As R increases →

V_out decreases faster →

takes less current for same brightness →

a more digital, less analog behavior

The extremes (R = 0 or R = inf) are tricky because they vary with op-amp characteristics. Theoretically, light always on with wire (0), always off with air (inf). Whether shine light or not!



recall: voltage

difference lights LED

Where are we/going?

Discuss goal of each phase/experiment with a partner

Construction phase:

Make Cph8 library (NK) and put into cells (109)

Prepare to identify mutants, TS41 X in Cylo
in the BP system, that might improve contrast

Evaluation phase:

DNA Sequencing what is the mutant? * other mystations?

(β-gal assay does the mutation seem to affect relative

protein

(β-gal) itself changed * * pokential for at hours *

(Western is Econol itself changed * * pokential for at hours *

(also: site distance of the mutant of the below

Other -- 12

Chance for design/investigation

- Day 7 = Western = lots of long incubations
- Chance to design and carry out an additional experiment that could shed light on your results
- Start thinking about it now (or rather as soon as you have your first batch of mutant data)
- Also a day to share, collaborate with other groups

Today in Lab: M2D5

- Miniprep each mutant candidate (light or dark)
- Resuspend in sterile water, mix an aliquot with sequencing primer → ideally done by 3 pm!
- β-gal assay: WT, HK dead, and your mutants
 - choose one dilution that should give credible data
 - perform in triplicate, towards statistics
- Prepare culture for protein analysis
- Leslie will discuss elements of a successful research article (esp. Discussion) @ 4 pm