

M2D1: Complete *in silico* cloning of dCas9 & actual confirmation digest

10/13/16

1. Pre-lab Discussion
2. Design primers to dCas9
3. *In silico* PCR amplification, digest, and ligation
4. *Actual* (wet-lab) diagnostic digest of pdCas9

(Almost) done with Mod1!

- Mini-presentation
 - due 10pm on Saturday, October 14
 - You don't have to be exhaustive;
tell a focused story
- Data summary
 - ✓ draft due 10pm on Wednesday, October 11
 - receive all comments by Monday, October 16th
 - revision due 10pm on Sunday, October 22nd
- Blog post
 - due 10pm on Monday, October 23



Sign up for journal club (M2D2 homework)

- Sign up on wiki for which day you will present:
M2D4 (October 25th) or M2D6 (November 1st)
- Pick 1 of 20 papers, or suggest your own
- Reserve paper by adding name next to it [Bagnall/WF/TeamColor]
 - First come first serve!
 - Only one T/R and one W/F student per article
 - Don't pick a paper randomly

Slot	Day 4 (T/R)	Day 6 (T/R)	Day 4 (W/F)	Day 6 (W/F)
1				
2				
3				
4				

“insert”
“fragment”
dCas9

“vector”
“backbone”
p_{tet}O-1

How is DNA engineered?

1. PCR amplification of DNA:

specific primers to gene of interest

dNTPs

Polymerase

2. Digestion:

Restriction enzymes (endonuclease)


Sticky (or blunt) ends

3. Ligation:

DNA Ligase

Seals backbone


Phosphodiester bonds



1. PCR amplification of DNA

top strand: 5'-->3'



Defining terminology



1. PCR amplification of DNA

Designing primers


- **Landing sequence:** match to dCas9
- **Flap sequence:** contains endonuclease recognition sequence and junk DNA



1. PCR amplification of DNA

Primer design guidelines

- Length
 - 17-28 base pairs
 - long enough to be specific, short enough for easy annealing
- GC content
 - 40-60%
 - AT: 2 H bonds
 - GC: 3 H bonds
 - GC clamp at ends
- T_m (primer)
 - $< 65^\circ\text{C}$
- Avoid secondary structures
 - hairpins
 - complementation w/in primer sequence
- Avoid repetitive sequences
 - Max of 4 di-nucleotide repeats (ex. ATATAT)
 - Max of 4 bp in a run (ex. GATGGGG)



1. PCR amplification of DNA

Three major PCR steps—which temperature & why?

- Melt
 - 95 °C
 - break hydrogen bonds
- Anneal depends on primer sequence
 - $T_m(\text{primer}) = 1/2 \text{ primer annealed to target}$
 - $T_{\text{anneal}} \sim T_m(\text{primer}) - 5^\circ\text{C}$
- Extend
 - 72 °C (for Taq)
 - 1000 bases/min

Leslie's favorite PCR animation

<http://learn.genetics.utah.edu/content/labs/pcr/>

cycle # 2

50° C

The temperature is lowered so the primers will attach.

<< BACK


NEXT >>

PCR

© 2008 GENETIC SCIENCE LEARNING CENTER, UNIVERSITY OF UTAH


2. Digestion

Restriction endonucleases create sticky ends on dCas9 insert and plasmid backbone





3. Ligation

Insert dCas9 into expression vector (backbone) to create new plasmid (pdCas9)



Confirmation digest considerations

- Do you have access to the enzymes?
- Are the two enzymes compatible?
- Are fragments easily distinguished on an agarose gel?



Today in lab

1. Reproduce *in silico* (in Benchling) the cloning of pdCas9
 - Design primers that would amplify the gene dCas9
 - Depict PCR amplification product
 - Digestion of dCAS9 PCR product and vector by restriction enzymes
 - Ligation of insert and vector ** at 2:30pm we will all work through this calculation together
2. Set up confirmation digests of pdCas9 for agarose gel electrophoresis
 - Choose restriction enzymes for diagnostic digest
 - Calculate volumes of digest components
 - Set-up digest and leave overnight at 37°C

~~at 2:30pm we will all work through this calculation together~~ ? For now, skip Part 3: #1-6

Ligation calculation

Goal: Calculate volumes of insert and backbone needed for ligation



Knowns:

- Need 50-100 ng backbone
- Backbone: 2592 bp
- Insert: 4113 bp
- Molar mass $\sim 660\text{g}/(\text{mol} \cdot \text{bp})$
- Desired molar ratio of insert to backbone is 4:1

Missing information needed to know what volumes of backbone and insert to use:

volume concentration of DNA

Use recovery gel to estimate insert and backbone concentrations



- Amount of backbone = 200 ng

- Amount of insert = 100 ng

- Concentration of backbone = $200 \text{ ng} / 5 \mu\text{L} = 40 \text{ ng}/\mu\text{L}$

- Concentration of insert = $100 \text{ ng} / 5 \mu\text{L} = 20 \text{ ng}/\mu\text{L}$

mass of DNA \neq molar amount of DNA

Goal: Calculate volumes of insert and backbone needed for ligation

Knowns:

- Need 50-100 ng backbone
- Backbone: 2592 bp
- Insert: 4113 bp
- Molar mass $\sim 660\text{g}/(\text{mol} \cdot \text{bp})$
- Desired molar ratio of insert to backbone is 4:1
- Concentration of insert: 20 ng/uL
- Concentration of backbone: 40 ng/uL

Calculate 4:1 (insert:backbone) *molar* amounts final volumes for ligation

1. Calculate moles of backbone

- $2592 \text{ bp} * (660 \text{ g} / (\text{mol} * \text{bp})) = 1.71 \times 10^6 \text{ g/mol}$
- so $50 \text{ ng} / (1.71 \times 10^6 \text{ g/mol}) = 2.9 \times 10^{-14} \text{ mol}$

2. Determine moles of insert needed (4x backbone)

- $4 \times 2.9 \times 10^{-14} \sim 1.2 \times 10^{-13} \text{ mol}$
- with $4113 \text{ bp} * (660 \text{ g} / (\text{mol} * \text{bp})) = 2.7 \times 10^6 \text{ g/mol}$
- so use $1.15 \times 10^{-13} \text{ mol} * 2.7 \times 10^6 \text{ g/mol} \sim 310 \text{ ng}$

1.2

3. Calculate volume of backbone and insert needed

- Backbone: $50 \text{ ng} / (40 \text{ ng/uL}) = 1.25 \text{ uL}$
- Insert: $310 \text{ ng} / (20 \text{ ng/μL}) = 15.5 \text{ μL}$

scale down: 1 uL backbone, 12.4 uL insert

Optimal backbone-to-insert ratio

- Ideally, want 4:1 insert:backbone
 - molar ratio, **not** mass or volume
- What if too much insert?
multiple inserts daisy-chained, inserted into backbone
- What if too much backbone?
backbone ligating to other backbone